DiaMetra

THYROGLOBULIN

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6. Pacini F., Elisei R., Fugazzola L., Pinchera A. Humoral Markers for Thyroid Carcinoma in Clinical Practice. *Diagn. Oncol.*, 1991, **1**, 194-196.
7. Chatherine Massart, Didier Maugendre. Importance of detection Method for Thyroglobulin Antibodies for the validity of Thiroglobulin measurements in sera from patients with Graves Disease. *Clin. Chem.* 2002, **48** (1), 108-107.

Ed 06/2010

DCM048-5

DiaMetra S.r.l. Headquater: Via Garibaldi, 18 – 20090 SEGRATE (MI)

Tel. 0039-02-2139184 – 02-26921595

Fax 0039-02-2133354.

Manufact: Via Giustozzi, 35/35a (già via Bartolomei) – Z.I Paciana – 06034 FOLIGNO (PG) ITALY.

Tel. 0039-0742-24851

Fax 0039-0742-316197

E-mail: info@diametra.com

DiaMetra

THYROGLOBULIN

for routine analysis

Determination of Thyroglobulin in human serum or plasma.

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO048

INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of Thyroglobulin concentration in human serum or plasma.

CLINICAL APPLICATIONS

Thyroglobulin (TG), a glycoprotein with a molecular weight of about 660,000 Daltons, is the thyroid's main iodine protein and the most important compound of follicular colloid. Thyroglobulin is the form under which the active hormones T₃ and T₄ together with their immediate forerunners MIT and DIT are laid inside the thyroid gland. The clinical applications of the TG dosage seem to originate from its specificity for the thyroid and related cells.

The dosage of TG can be used as support to scintigraphies or other techniques for studying pathogenesis, making a diagnosis and analyzing the course of thyroid disorders.

The dosage of TG before and after replacement treatment with L-Thyroxin cannot be established in cases of hypothyroidism due to thyroid agenesis. In cases of secondary hypothyroidism with a dysglandular goiter or ectopic thyroid, the levels of TG are normal or high. The circulating levels of TG tend to increase in several thyroid disorders such as toxic and atoxic goiter, subacute thyroiditis, Basedow's disease and carcinoma. In Basedow's disease the TG dosage is a potentially interesting index of normalization of hyperthyroidism in patients treated with anti-thyroid drugs. In the oncology field and more specifically for differentiated thyroid carcinoma, there are very promising applications linked to the ability of thyroid tumors tissues to concentrate iodine and synthesize TG as a normal thyroid. Basically the dosage of TG can be used as follows:

a. Pre-operating diagnosis of thyroid tumors.

This application does not allow the differentiated diagnosis of the tumor as the values of TG seen in malignant and benign nodules are superimposable.

b. Post-operation monitoring

In patients treated surgically or with radiotherapy, long lasting TG levels suggest the presence of a residual carcinoma and/or carcinoma with metastasis.

c. Monitoring of totally thyroidectomized patients

The use of circulating TG as an indicator of recurrent tumors (metastasis marker) has an established clinical value: the increase of Thyroglobulinaemia indicates the need to undergo further analysis for confirming the diagnosis. Interesting advantages can come from: a) a reduced use of scintigraphic diagnostic techniques as they imply regular suspension of replacement treatment and frequent exposure to radiation, b) and complete completion of the information obtained via scintigraphy

1. PRINCIPLE

This kit is based on an immunoenzymometric assay (IEMA). Four different anti-TG monoclonal antibodies are used, three coated on the wells and the fourth conjugated to Biotin. During incubation, the TG present in the calibrators and samples is bound to the antibodies and biotinylate forming a "sandwich". At the end of the incubation, the unbound material is removed by an aspiration/washing cycle. A solution of streptavidin conjugated with peroxidase (HRP) is then added to all wells. This will react with the biotinylated antibodies which are bound to the well. After a further aspiration/washing cycle, the residual enzyme activity found in the wells will thus be directly proportional to TG concentration in the calibrators and samples and evidenced by adding into the wells a Chromogen solution (tetramethylbenzidine, TMB) in a Substrate-Buffer. Colorimetric reading will be performed by using a spectrophotometer at 450 nm and at 405 nm wavelength.

2. REAGENTS, MATERIALS AND INSTRUMENTATIONS

2.1 Reagents and materials supplied in the kit

1. Reference Serum- (Standard)

STD ₀ (1 vial) 2.0 mL	REF DCE002/4806-0
STD ₁ (1 vial) 1.0 mL	REF DCE002/4807-0
STD ₂ (1 vial) 1.0 mL	REF DCE002/4808-0
STD ₃ (1 vial) 1.0 mL	REF DCE002/4809-0
STD ₄ (1 vial) 1.0 mL	REF DCE002/4810-0
STD ₅ (1 vial) 1.0 mL	REF DCE002/4811-0
STD ₆ (1 vial) 1.0 mL	REF DCE002/4812-0
	REF DCE021/4821-0

2. Control serum I

TG in proteic matrix (1 vial) 1 mL

3. Control serum II

TG in proteic matrix (1 vial) 1 mL

4. Anti humanTg Biotin

Biotinylated anti human TG (1 bottle) 13 mL

5. Enzyme Conjugate (1 bottle) 15 mL

Streptavidine-HRP

6. Coated Microplate

Anti TG coated microplate (1 micropatle breakable)

7. Recovery solution (1 vial) 3 mL

TG in proteic matrix 50 ng/mL

8. Concentrated wash solution 20X (1 bottle) 50 mL

NaCl 9 g/L; Tween-20 22 g/L

9. TMB-substrate (1 bottle) 12 mL

H₂O₂-TMB 0.26 g/L (avoid any skin contact)

10. Stop solution (1 bottle) 12 mL

Sulphuric acid 0.15 mol/L (avoid any skin contact)

2.2 Note

- Store the reagents at 2±8°C
- Open the bag of reagent 6 (Coated Microplate) only when it is at room temperature and close immediately after use. Once opened, it is stable until the expiration date of the kit
- Do not use reagents beyond the kit expiration date.
- Opened reagents are stable for sixty (60) days when stored at 2±8°C.
- Human Serum Reference (STANDARD) 7 vials of serum reference for at concentrations of 0, 1, 3, 10, 30, 100, 300 ng/mL. Store at 2±8°C. A preservative has been added.

MATERIALS Required But Not Provided:

1. Pipette capable of delivering 50µl volumes with a precision of better than 1.5%.
2. Dispenser(s) for repetitive deliveries of 0.100mL and 0.300mL volumes with a precision of better than 1.5%.
3. Microplate washer or a squeeze bottle (optional).
4. Microplate Reader with 405 nm, 450nm and 620nm wavelength absorbance capability.
5. Absorbent Paper for blotting the microplate wells.
6. Vacuum aspirator (optional) for wash steps.
7. Timer.

3.WARNINGS AND PRECAUTIONS

In order to obtain correct and reproducible results, the following rules must be observed:

- Do not mix specific reagents from different lots.
- It is possible to use common reagents from different lots.
- Do not use reagents beyond their expiry date.
- Do not store or leave reagents and samples at high temperatures or areas of possible contamination.
- Use thoroughly clean glassware, free from metal ion contamination or oxidizing substances.
- Use distilled or deionized water, stored in perfectly clean containers.
- Carefully avoid any contamination among samples; for this purpose, disposable tips should be used for each sample and reagent.
- Do not modify in any way the "Assay Procedure". If you not respect:
 - exact incubation times and quantities adding the reagents incubation times and temperature may cause incorrect clinical results.
 - Reconstitute lyophilized reagents, if present, as described on the relative labels. Any deviation in reagent use or wrong volumes, may affect the reliability of results obtained.
 - In case of manual procedure, it is important to use calibrated pipettes and have appropriate technical manuals. Primary importance is a good precision preparing and dispensing the reagents. Ensure that all the equipment used is in perfect working order, has been correctly calibrated and is regularly maintained.
 - Ensure that all the equipment used (glassware, dry heater, microplate shaker, microplate washers, spectrophotometer and fridge/freezers used for reagent and sample storage) is in perfect working order, has been correctly calibrated and is regularly maintained. Any deviation from the correct use of the equipment listed can produce errors in the methodology, this may affect the reproducibility and reliability of results obtained.

- Utilise a suitable method for the correct identification of patient samples. Incorrect identification may cause a specificity losses of the system and wrong clinical results.

In order to avoid personal and environmental contamination, the following precautions must be observed:

- Use disposable gloves while handling potentially infectious material and while performing the assay.
- Do not pipette reagents by mouth.
- Do not smoke, eat, drink or apply cosmetics during the assay.
- Chromogen and Blocking Reagent should be handled with care. Avoid contact with skin, eyes and mucous membranes. In case of accident rinse thoroughly with running water.
- All material of human origin used for the preparation of this kit tested negative for HBsAg, anti-HIV and anti-HCV. Since no test at present can guarantee complete absence of these viruses, all samples and reagents containing biological material used for the assay must be considered potentially infectious.
- Avoid splashing and aerosol formation; in such cases, carefully wash with a 3% sodium hypochlorite solution. Any such cleaning material must be treated as potentially infectious and disposed of accordingly.
- Some reagents contain sodium azide as preservative; to prevent build-up of explosive metal azides in lead and copper plumbing, reagents should be discarded by flushing the drain with large amounts of water.
- The reagents for which a Safety Data Sheet is not supplied do not contain hazardous chemical substances or, if they do, these are below the concentration limits established by the Italian decree D.Lgs.285/98 in compliance with EEC directive 91/155.
- According to Italian decree D.L. no. 22 dated 05.02.97, in compliance with EEC directives (91/156/EEC, 91/689/EEC, 94/62/EEC), all waste products originating from either manual and/or automated processing are classified as hazardous special waste material (European classification code180103). As such, they must be eliminated by delegating to special enterprises, qualified for waste collection and disposal.

3.1 REAGENT PREPARATION:

Wash Solution - Dilute the contents of each vial of the buffered wash solution concentrate (20x) with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:20 dilution ratio. The diluted wash solution is stable for 30 days at 2±8°C.

4. SPECIMEN COLLECTION AND PREPARATION

The assay can be performed in human serum or plasma samples. Highly lipemic or hemolyzed samples must be discarded. Keep samples at 2±8°C for 1-2 days; for longer periods it is advisable to freeze samples at -20°C. Repeated freezing and thawing of samples should be avoided. Samples with TG concentrations higher than 300 ng/mL, must be diluted with the Zero Calibrator. We suggest a dilution of 1:5 (100 µL of sample + 400 µL of Zero Calibrator).

5 PROCEDURE

As it is necessary to perform in duplicate, each determination should also include two wells for each *Sample*, two wells for each *Serum Control* and two wells and for each point of standard curve.

	Blank	Standard	Sample/Control
Sample/Control	---	---	50 µL
Standard S ₀ – S ₆	---	50 µL	---
Anti Tg biot.	---	100 µL	100 µL

Incubate at 37°C for 90 minutes

Remove the contents from each well.

Wash three times each well with 0.3 mL of diluted washing solution, leave the excess liquid to drain away by inverting the plate on absorbent paper. Dispense into each well:

Dispense:

	Blank	Standard	Sample/Control
Enzyme conjugate	---	100 µL	100 µL

Incubate at 37°C for 30 minutes

Remove the contents from each well.

Wash three times each well with 0.3 mL of diluted washing solution, leave the excess liquid to drain away by inverting the plate on absorbent paper. Dispense into each well:

	Blank	Standard	Sample/Control
TMB Substrate	100 µL	100 µL	100 µL

Incubate at Room Temperature (22-28°C) for 15 minutes in the dark.

Dispense

	Blank	Standard	Sample/Control
Stop solution	100 µL	100 µL	100 µL

Shake the microplate gently. Determine the absorbance (E) by reading at 450 nm for each test-well (within 30 minutes after the addition of the Stop solution) after zeroing the instrument with the S₀.

6. CALCULATION OF RESULTS

In order to obtain a better sensitivity, the present method employs spectrophotometric reading at two wavelengths (450 and 405 nm). For samples with TG concentrations ranging from 0 to 30 ng/mL, read at 450 nm wavelength; for samples with TG level higher than 30 ng/mL, read at 405 nm wavelength.

Draw a calibration curve on millimetric graph paper, by plotting the calibrator concentrations (x-axis) against the absorbance obtained for each calibrator (y-axis). Corresponding TG concentrations in µIU/mL are obtained by interpolating the absorbances of each sample on the calibration curve; in case of diluted samples, multiply by the dilution factor.

6.1 Example of Calculation

The values shown below must be considered as an example and must not be used in place of experimental data.

Standard/Sample	O.D. 450 nm	TG	O.D. 405 nm	TG
Standard 0 ng/mL	0.020		0.007	
Standard 1 ng/mL	0.047		0.015	
Standard 3 ng/mL	0.127		0.042	
Standard 10 ng/mL	0.313		0.098	
Standard 30 ng/mL	0.863		0.373	
Standard 100 ng/mL	1.650		0.603	
Standard 300 ng/mL	> 3.000		1.488	
Sample 1	0.450	15.4 ng/mL	0.141	
Sample 2	2.240		0.829	166.4 ng/mL

6.2 Normal Values

The values reported below are indicative. We suggest that each laboratory establishes its own normal range.

Normal values have been established on 85 healthy subjects without thyroid disorders. Thyroglobulin levels were **lower than 40 ng/mL**.

6.3 Validation Criteria

Before proceeding to calculation of results, make sure the control serum concentration is included within the value described on the Quality Control Sheet.

7. PERFORMANCES OF THE ASSAY

7.1 Specificity

No cross reaction was observed with MIT, DIT, rT3, T3, T4, TSH, FSH and LH. This analytic method showed a cross-reactivity of 0.01% with human TBG.

7.2 Sensitivity

The sensitivity was calculated upon the calibration curve and expressed as the minimal dose showing a significant difference from the Zero Calibrator (mean value + 2 S.D.). This dose is **0.15 ng/mL**.

7.3 Precision

Precision was evaluated determining the repeatability and the reproducibility of the assay (intra- and inter-assay variability), on 3 sera at different TG concentrations.

Repeatability (Intra-assay)

Serum	Media (ng/mL)	±	D.S.	% C.V.	Replicates n.
1	152.41	±	5.71	3.75	10
2	15.14	±	0.26	1.71	10
3	1.69	±	0.10	6.14	10

Reproducibility (Inter-assay)

Serum	Media (ng/mL)	±	D.S.	% C.V.	Replicates n.
1	150.4	±	9.3	6.2	10
2	38.1	±	3.0	7.8	10
3	1.8	±	0.13	6.9	10

7.4 Accuracy

Accuracy of the method has been checked by the recovery and parallelism tests:

Recovery Test.

Known amounts of TG have been added to two normal sera and tested.

Added (ng/mL)	Measured (ng/mL)	Recovered (ng/mL)	Recovery %
S1	1.9	---	---
S1 + 12.5	14.9	13.0	104.0
S1 + 25	25.2	23.3	93.2
S1 + 50	50.7	48.8	97.6
S1 + 100	93.6	91.7	91.7
S1 + 200	195.00	193.1	96.6
S2	17.4	---	---
S2 + 12.5	27.0	9.6	76.8
S2 + 25	42.0	24.6	98.4
S2 + 50	65.0	47.6	95.2
S2 + 100	116.2	98.8	98.8
S2 + 200	205.1	187.7	93.9

Parallelism Test.

Two sera with high TG concentration were tested at different dilutions with the Zero Calibrator.

Dilution	Expected (ng/mL)	Measured (ng/mL)
S1 undiluted	---	179.50
1:2	89.75	84.75
1:4	44.88	47.57
1:8	22.44	28.46
S2 undiluted	---	25.69
1:2	12.85	15.05
1:4	6.42	6.97
1:8	3.21	3.83

7.5 Correlation

The Diametra TG ELISA Kit was compared to another commercially available TG assay (TG Zentech Irma Kit). 66 serum samples were analysed according in both test system. The linear regression curve is:

$$TG \text{ Zentech} = 1.029 * TG \text{ Diametra} + 2.069 \quad (R^2=0.952)$$

8. RECOVERY INTO THE SERUM SAMPLE

In a sample, the presence of anti-TG self-antibodies influences TG assay and thus can cause inaccurate results. Therefore, it is necessary to carry out a clinical recovery test to confirm the accuracy of the result. This test should not be considered as a method for discovering anti-TG antibodies.

Procedure:

Dilute 1/2 of the serum sample with the Recovery solution, e.g. 50 µL of sample + 50 µL of recovery solution. The TG concentration in the recovery solution is 50 ng/mL (SR).

Test the undiluted sample (S1) and the sample half diluted with the recovery solution according to the Assay scheme.

The percentage of recovered TG for a sample is calculated as follows:

$$\text{Recovery (\%)} = \frac{\text{ng/mL sample S2}}{(\text{ng/mL sample S1} + 50)/2} \times 100$$

Recoveries less than 75% or more than 120% indicate interferences from anti-TG self-antibodies.

12. LIMITS OF THE ASSAY

The results of the assay must be carefully interpreted and confirmed by clinical evaluations and further diagnostic tests.

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
ERROR POSSIBLE CAUSES / SUGGESTIONS

No colorimetric reaction

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation