DiaMetra

TSH RECEPTOR AUTOANTIBODY

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BIBLIOGRAFIA

- J. Bolton et al Measurement of thyroid stimulating hormone receptor autoantibodies by ELISA Clin. Chem 1999 45: 2285-2287
- K Kamijo TSH receptor antibody measurement in patients with various thyrotoxicosis and Hashimoto's thyroiditis: a comparison of two two-step assays, coated plate ELISA using porcine TSH receptor and coated tube radioassay using human recombinant TSH receptor Endocrine Journal 2003 50:113-116
- B. Rees Smith et al A new assay for thyrotropin receptor autoantibodies Thyroid 2004 14: 830-835

Ed. 02/2011

DCM085-5

DiaMetra S.r.l. Headquarter: Via Garibaldi, 18
20090 SEGRATE (MI) Italy
Tel. 0039-02-2139184 – 02-26921595
Fax 0039-02-2133354.

Manufactory: Via Giustozzi 35/35a – Z.I Paciana
06034 FOLIGNO (PG) Italy
Tel. 0039-0742-24851
Fax 0039-0742-316197
E-mail: info@diametra.com

DiaMetra

TSH RECEPTOR AUTOANTIBODY

for routine analysis

Quantitative determination of thyrotropin receptor autoantibodies (TRAb) in human serum

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO085

INTENDED USE

TSH receptor autoantibody (TRAb) ELISA kit is intended for use by professional persons only for the quantitative determination of thyrotropin receptor autoantibodies in human serum.

TRAb ELISA kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

The hyperthyroidism is a disease that produce an increase in blood circulating levels of thyroid hormones. In some cases, as in Graves' disease is caused by autoantibodies against the TSH receptor (TSHr). Autoantibodies mimic the effect of TSH on the thyroid gland causing a rise in blood levels of T3 and T4. The determination of these autoantibodies (TRAb) may be useful for the diagnosis of the disease and its treatment. The main treatment for Graves' disease is represented by use of antithyroid drugs (propylthiouracil or methimazole), I^{131} and surgery. The dosage of TRAb is therefore useful in the course or at the end of therapy.

2. PRINCIPLE

In TRAb ELISA TSH receptor autoantibodies patient sera, calibrators and controls are allowed to interact with TSH receptor coated into ELISA plate wells. After a 2 hours incubation, the samples are discarded leaving TRAb bound to the immobilized TSH receptor. TSH Biotin, added in a 2nd incubation step, interacts with immobilized TSH receptors, which have not been blocked by the bound TRAb from patient sera, calibrators or controls. The amount of TSH Biotin bound to the plate is then determined in a third incubation step by addition of streptavidin peroxidase, which binds specifically to biotin. Excess unbound streptavidin peroxidase is then washed away and the addition of tetramethylbenzidine (TMB) substrate results in formation of a blue colour.

This reaction is stopped by the addition of stop solution causing the well contents to turn from blue to yellow. The absorbance of the yellow reaction mixture at 450nm is then read using an ELISA plate reader. A lower absorbance indicates the presence of TRAb in the test sample as TRAb inhibits the binding of TSH biotin to TSH receptor coated plate wells.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

- TSH Receptor Standards (4 vials, 1 mL each)
 - STD1 **REF** DCE002/8507-0
 - STD2 **REF** DCE002/8508-0
 - STD3 **REF** DCE002/8509-0
 - STD4 **REF** DCE002/8510-0
- Controls (2 vials, 1mL each, ready to use)
 - Negative Control **REF** DCE045/8501-0
 - Positive Control **REF** DCE045/8502-0
- TSH biotin (3 vials, lyophilized) **REF** DCE019/8519-0
- TSH biotin buffer (1 vial, 15 mL) **REF** DCE047/8547-0
- 20X Conc. Streptavidin-Peroxidase (1 vial, 0.75 mL) **REF** DCE041/8541-0
- Streptavidin-Peroxidase Diluent (1 vial) 15 mL **REF** DCE048/8548-0
- Coated Microplate
(1 breakable microplate coated with TSH receptor) **REF** DCE002/8503-0
- Start Buffer (1 vial, 10 mL) **REF** DCE046-0
- TMB Substrate (1 vial, 15 mL) **REF** DCE004/8504-0
- 10X Conc. Wash Solution (1 vial, 100 mL) **REF** DCE006/8506-0
- Stop Solution (1 vial, 10 mL)
0.25 M sulfuric acid **REF** DCE005/8505-0

3.2. Reagents necessary not supplied

Distilled or deionized water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm)

Note

After opening, return any unused wells to the original foil packet and seal. Then place foil bag in the self-seal plastic bag with desiccant provided. Store at 2-8°C for up to 6 months.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only.
- Use appropriate personal protective equipment while working with the reagents provided.
- All human source material used in the preparation of standards and controls for this product has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Standard and the Controls should be handled in the same manner as potentially infectious material.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of Sodium Azide (NaN₃) or Proclin 300^R as preservatives. Avoid the contact with skin or mucosa.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, allow scroll through large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry dates printed on the labels of the box and of the vials must be observed. Do not use any kit component beyond their expiry date.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is

advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinization and decontamination are effective in avoiding contamination of the conjugate; this procedure is highly recommended when the kit is processed using analyzers which are not equipped with disposable tips.

For this purpose, Diametra supplies a separate decontamination reagent for cleaning needles.

- If you use automated equipment is your responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated should not be used in the assay. Highly lipemic or haemolysed specimens should similarly not be used
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Standard (S₁,S₂,S₃,S₄) Concentration

The standard calibrated against the WHO NIBSC 90/672 have approximatively the following concentration:

	S1	S2	S3	S4
U/L	1	2	8	40

6.2. Preparation of the Sample

Sera to be analysed should be assayed soon after separation or stored, preferably in aliquots, at -20°C or below.

Repeated freeze thawing or increases in storage temperature must be avoided.

Incorrect storage of serum samples can lead to loss of TRAb activity. Do not use highly lipaemic or haemolysed serum samples. Do not use plasma in the assay. When required, thaw test sera at room temperature and mix gently to ensure homogeneity.

The Controls are ready to use.

6.3. Preparation of the TSH biotin

Reconstitute each vial with 4.5 mL TSH Biotin buffer. When more than one vial is to be used, pool the vials and mix gently before use. Store at 2-8°C for up to 6 months after reconstitution.

6.4. Preparation of the Streptavidin Peroxydase

Dilute 1:20 with Streptavidin Peroxydase Diluent (for example, 0.5 mL of Streptavidin-Peroxidase 20X concentrate + 9.5 mL of Streptavidin-Peroxidase Diluent). Store at 2-8°C for up to expiry date of the kit.

6.5. Preparation of the Wash Solution

Prepare a sufficient amount of washing solution by diluting the concentrated wash buffer 1:10 with distilled or de-ionized water. For example, dilute 50 mL of the concentrate with 450 mL of distilled water. The solution should be free of crystals before dilution, otherwise dissolve by warming up to max. 37°C. The diluted washing solution can be stored at 2-8 °C up to expiry date of the kit.

6.6. Procedure

Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes. As it is necessary to perform the determination in duplicate, prepare two wells for each of the four points of the standard curve (S₁-S₄), two for each sample, one for Blank.

TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22-28°C) for 30 minutes in the dark.			
Stop Solution	50 µL	50 µL	50 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against Blank within 20 min after adding the stop solution.			

7. QUALITY CONTROL

It is recommended that each laboratory include its own panel of control samples in the assay. Each laboratory should establish its own normal and pathological reference ranges for TRAb levels.

8. RESULTS

8.1. Standard Curve

A calibration curve can be established by plotting calibrator concentration on the x-axis (log scale) against the absorbance of the calibrators on the y-axis (linear scale). The TRAb concentrations in patient sera can then be read off the calibration curve. Other data reduction systems can be used. Results can also be expressed as inhibition (%) of TSH binding calculated using the formula;


$$100 \times \left(1 - \frac{\text{test sample absorbance at 450 nm}}{\text{negative control (D1) absorbance 450 nm}} \right)$$

Samples with high TRAb concentrations can be diluted in kit negative control (D1). For example, 20 µL of sample plus 180 µL of negative control to give a 10x dilution. Other dilutions (e.g. 100x) can be prepared from a 10x dilution or otherwise as appropriate. Some sera will not dilute in a linear way and we suggest that the dilution giving a value closest to 50% inhibition is used for calculation of TRAb concentration.

Typical Results (example only, not for use in calculation of actual results)

Sample	A450	% Inhibition	U/L
Control D1	2.00	0	0
C1	1.70	15	1
C2	1.50	25	2
C3	0.65	68	8
C4	0.15	93	40
Control D2	1.26	37	3.5

Reagents	Standard	Sample or Controls	Blank
Start buffer	75 µL	75 µL	
Standard	75 µL		
Sample or Controls		75 µL	
Cover the plate and incubate 120 minutes at room temperature (22-28°C) while shaking at >500 rpm. Alternatively incubate without shaking for at least 180 minutes at room temperature (22-28°C). Remove the contents from each well and wash the wells with 300 µL of diluted Wash Solution, drain the wash completely.			
TSH Biotin	100 µL	100 µL	
Cover the plate and incubate 25 minutes at room temperature (22-28°C). Remove the contents from each well and wash the wells with 300 µL of diluted Wash Solution, drain the wash completely.			
Streptavidin-Peroxidase	100 µL	100 µL	
Cover the plate and incubate 20 minutes at room temperature (22-28°C). Remove the contents from each well and wash the wells with 300 µL of diluted Wash Solution. Repeat the washing procedure other two times by draining the wash completely.			

8.2. Reference values

anti-TRAb	
negative	≤ 1.0 U/L
Grey zone	1.1 – 1.5 U/L
positive	> 1.5 U/L

NB: The above mentioned reference values provide only a guide.

9. PERFORMANCE AND CHARACTERISTICS

9.1. Clinical Specificity

154 samples from healthy blood donors were assayed in the TRAb ELISA kit. 152 (99%) were identified as being negative for TSH Receptor autoantibodies.

9.2. Clinical Sensitivity

50 samples from patients diagnosed with Graves' disease were assayed using the TRAb ELISA kit. 49 (98%) were identified as being positive for TSH Receptor autoantibodies. 1 sample (2%) was identified as being within the equivocal range.

9.3. Lower Detection Limit

The kit negative control was assayed 32 times and the mean and standard deviation calculated. The lower detection limit at 2 standard deviations was 0.21 U/mL.

9.4. Intra and inter-assay variations

9.4.1 Intra-Assay

Within run variation was determined by replicate 16 times two different sera with values in the range of standard curve. The within assay variability is ≤ 7.6%

9.4.2 Inter-Assay

Between run variation was determined by replicate the measurements of one control serum with different lots of kits and/or different mix of lots of reagents. The between assay variability is ≤ 6.5%.

9.5. Clinical Accuracy

Analysis of sera from patients with autoimmune diseases other than Graves' disease indicated non interference from autoantibodies to: thyroglobulin, thyroid peroxidase, glutamic acid decarboxylase, 21-hydroxylase, acetylcholine receptor, dsDNA or from rheumatoid factor.

9.6. Interference

No interference was observed when samples were spiked with the following materials: haemoglobin (5 mg/mL), bilirubin (20 mg/dL), intralipid (10 mg/mL).

10. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

BIBLIOGRAPHY

- J. Bolton et al Measurement of thyroid stimulating hormone receptor autoantibodies by ELISA Clin. Chem 1999 45: 2285-2287
- K Kamijo TSH receptor antibody measurement in patients with various thyrotoxicosis and Hashimoto's thyroiditis: a comparison of two two-step assays, coated plate ELISA using porcine TSH receptor and coated tube radioassay using human recombinant TSH receptor Endocrine Journal 2003 50:113-116
- B. Rees Smith et al A new assay for thyrotropin receptor autoantibodies Thyroid 2004 14: 830-835

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ERROR POSSIBLE CAUSES / SUGGESTIONS

No colorimetric reaction

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation