DiaMetra

Anti dsDNA IgG

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BIBLIOGRAFIA

1. Condemi, Jhon J. The autoimmune disease. The Journal of the American Medical Association 1987, Vo 258 n 20 2920 - 2929
2. Reichlin, M., Van Venrooij, W. J.: Autoantibodies to UNRP particles: relationship to clinical diagnosis and nephritis. Clin. exp. Immunol. 1991; 83:286-90.
3. Tan EM, et al.: The 1982 Revised Criteria for the classification of Systemic Lupus Erythematosus. Arthritis and Rheumatism 25: 1271-1277, 1982.
4. McCarty GA. Autoantibodies and their relation to rheumatic diseases. Medical Clinics of North America 70: 237-261, 1986.
5. Hardin JA. The lupus autoantigens and the pathogenesis of systemic lupus erythematosus. Arth Rheum 29: 457-460, 1986.
6. The first international standard for antibody to double stranded DNA. Annals of the Rheumatic Disease 198; Vo 47: 740 - 746
7. Homburger HA, et al: detection of Antinuclear antibodies, comparative evaluation of enzyme immunoassay and indirect immunofluorescence methods. Arch Pathol Lab Med 122:993-999 (1998)
8. Smeenk, R. et al. Avidity of Antibodies to dsDNA. Comparison of IFT on Crithidia Luciliae, FARR assay and PEG assay. The journal of Immunology 1982 Vo 128 n.1. 73 -78

Ed. 02/2011

DCM095-4

DiaMetra S.r.l. Headquater: Via Garibaldi, 18 –
20090 SEGRATE (MI) Italy
Tel. 0039-02-2139184 – 02-26921595
Fax 0039-02-2133354.

Manufact: Via Giustozzi 35/35a, Z.I Paciana – 06034
FOLIGNO (PG) Italy
Tel. 0039-0742-24851
Fax 0039-0742-316197
E-mail:info@diametra.com

DiaMetra

Anti dsDNA IgG

for routine analysis

Quantitative determination of IgG class antibodies against dsDNA in human serum or plasma

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO095

INTENDED USE

Anti dsDNA IgG kit is an indirect solid phase enzyme immunometric assay (ELISA) designed for the quantitative measurement of IgG class antibodies directed against dsDNA in human serum or plasma.

Anti dsDNA IgG is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Anti-dsDNA IgG test is used for initial diagnosis of Systemic Lupus Erythematosus (SLE) and for diagnosis of SLE different diseases. Besides the determination of high titers of antinuclear antibody (ANA), the determination of autoantibodies against dsDNA is one of the ACR criteria (American College of Rheumatology) for the diagnosis of Systemic Lupus Erythematosus (SLE). The determination of the concentration of antibodies can be used to monitor treatment success and predict possible attacks of the disease (SLE)

2. PRINCIPLE

Anti-dsDNA IgG test is based on the binding of serum or plasma IgG antibodies on dsDNA coated on the microplates. The antibodies in calibrators, controls or prediluted patient samples bind into the inner surface of the wells. After a 30 minutes incubation the microplate is washed with wash buffer for removing non-reactive serum components.

An anti-human-IgG horseradish peroxidase conjugate solution recognizes the IgG class antibodies bound to the immobilized dsDNA antigens. After a 30 minutes incubation any excess of enzyme conjugate, which is not specifically bound, is washed away with wash buffer.

A chromogenic substrate solution containing TMB is dispensed into the wells. After 15 minutes of incubation the color development is stopped by adding the stop solution. The solution color changes into yellow. The amount of color is directly proportional to the concentration of the anti-dsDNA IgG antibodies present in the original sample.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

- Anti dsDNA Standards (5 vials, 1.2 mL each)
Phosphate buffer 0.1 M, NaN₃ < 0.1%, human serum

STD0	REF DCE002/9506-0
STD1	REF DCE002/9507-0
STD2	REF DCE002/9508-0
STD3	REF DCE002/9509-0
STD4	REF DCE002/9510-0

- Controls (2 vials, 1.2 mL each, ready to use)
Phosphate buffer 0.1 M, NaN₃ < 0.1%, human serum

Negative Control	REF DCE045/9501-0
Positive Control	REF DCE045/9502-0

- Sample Diluent (1 vial, 100 mL)
Phosphate buffer 0.1 M, NaN₃ < 0.1%
REF DCE053-0
- Conjugate (1 vial, 15 mL)
Anti h-IgG conjugated with peroxidase, BSA 0.1%, Proclin < 0.0015%
REF DCE002/9502-0
- Coated Microplate
(1 breakable microplate coated with dsDNA)
REF DCE002/9503-0
- TMB Substrate (1 vial, 15 mL)
3,3',5,5'-tetramethylbenzidine 0.26 g/L, hydrogen peroxide 0.05%
REF DCE004-0
- 10X Conc. Wash Solution (1 vial, 50 mL)
Phosphate buffer 0.2M, proclin < 0.0015%
REF DCE054-0
- Stop Solution (1 vial, 15 mL)
Sulphuric acid 0.15 M
REF DCE005-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.
Microplates reader (450 nm)

4. WARNINGS

- This kit is intended for in vitro use by professional persons only.
- Use appropriate personal protective equipment while working with the reagents provided.
- All human source material used in the preparation of standards and controls for this product has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Standard and the Controls should be handled in the same manner as potentially infectious material.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of Sodium Azide (NaN₃) or Proclin 300^R as preservatives. Avoid the contact with skin or mucosa.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents,

allow scroll through large amounts of water to prevent azide build-up.

- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry dates printed on the labels of the box and of the vials must be observed. Do not use any kit component beyond their expiry date.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinization and decontamination are effective in avoiding contamination of the conjugate; **this procedure is highly recommended when the kit is processed using analyzers which are not equipped with disposable tips.** For this purpose, Diametra supplies a separate decontamination reagent for cleaning needles.
- If you use automated equipment it is your responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.

- Samples microbiologically contaminated should not be used in the assay. Highly lipemic or haemolysed specimens should similarly not be used
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. STORAGE CONDITIONS

- Store all the kit reagents at 2-8°C. Do not freeze. Reagents are stable until the expiration date when stored and handled as directed.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.

7. PROCEDURE

7.1. Preparation of the Standard (S₀,...,S₄)

The assay system is calibrated in relative arbitrary units. The standards have approximatively the following concentration:

	S ₀	S ₁	S ₂	S ₃	S ₄
AU/mL	0	15	30	60	240

7.2. Preparation of the Sample

For determination of Anti-dsDNA antibodies, human serum or plasma are the preferred sample matrices.

All serum and plasma samples have to be prediluted with sample diluent 1 : 100. Therefore 10 µL of sample may be diluted with 990 µL of sample diluent.

The patients need not to be fasting, and no special preparations are necessary. Collect blood by venipuncture into vacutainers and separate serum (after clot formation) or plasma from the cells by centrifugation.

Samples may be stored refrigerated at 2 - 8 °C for at least 5 days. For longer storage of up to six months samples should be stored frozen at -20°C. To avoid repeated thawing and freezing the samples should be aliquoted.

Neither Bilirubin nor Hemolysis have significant effect on the procedure.

The Controls are ready to use.

7.3. Preparation of the Wash Solution

Dilute the contents of each vial of the buffered wash solution concentrate (10x) with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

7.4. Procedure

Allow all reagents to reach room temperature (22-28°C).

As it is necessary to perform the determination in duplicate, prepare two wells for each of the five points of the standard curve (S₀-S₄), two for each Control, two for each sample, one for Blank.

Reagent	Standard	Sample or Controls	Blank
Standard S ₀ -S ₄	100 µL		
Controls		100 µL	
Diluted Sample		100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells three times with 300 µL of diluted wash solution			
Conjugate	100 µL	100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells three times with 300 µL of diluted wash solution.			
TMB substrate	100 µL	100 µL	100 µL
Incubate 15 minutes in the dark at room temperature (22-28°C).			
Stop solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against Blank.			

8. RESULTS

8.1. Standard curve

For Anti-dsDNA IgG a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice. Smoothed-Spline Approximation and log-log coordinates are also suitable. However we recommend using a Lin-Log curve. First calculate the averaged optical densities for each calibrator well. Use lin-log graph paper and plot the averaged optical density of each calibrator versus the concentration. Draw the best fitting curve approximating the path of all calibrator points. The calibrator points may also be connected with straight line segments. The concentration of unknowns may then be estimated from the calibration curve by interpolation.

Typical Results (example only)

The figure below shows typical results for Anti-dsDNA IgG. These data are intended for illustration only and should not be used to calculate results from another run.

N	OD1	OD2	mean	C1	C2	mean	CV%
STD0	0.007	0.008	0.008	0.00	0.03	0.01	---
STD1	0.426	0.403	0.415	15.90	14.97	15.44	4.28
STD2	0.731	0.717	0.724	29.68	28.98	29.33	1.68
STD3	1.208	1.243	1.226	59.17	61.90	60.53	3.19
STD4	2.272	2.218	2.245	251.3	229.3	240.3	6.48

9. REFERENCE VALUES

In a normal range study with serum samples from healthy blood donors the following ranges have been established with the Anti-dsDNA tests:

Antibodies anti-dsDNA IgG	[AU/mL]
normal:	< 30
elevated:	> 30

Positive results should be verified concerning the entire clinical status of the patient. Also every decision for therapy should be taken individually.

It is recommended that each laboratory establishes its own normal and pathological ranges of serum Anti-dsDNA.

9.1 Specificity

Comparison test against a commercial reference kit, performed on 36 sera (15 of them positive sera and 21 negative sera) showed a 95.24% specificity.

9.2 Sensitivity

Comparison test against a commercial reference kit, performed on 36 sera (15 of them positive sera and 21 negative sera) showed a 93.33% sensibility.

9.3 Detection limit

The lowest concentration of anti-dsDNA IgG that can be distinguished from zero standard is 0.42 AU/mL with a confidence limit of 95%.

9.4 Precision and Reproducibility

9.4.1 Intra-Assay

Within run variation was determined by replicate 16 times two different sera with values in the range of standard curve. The within assay variability is ≤ 3.4%

9.4.2 Inter-Assay

Between run variation was determined by replicate measurements of two different control sera with different lots of kits and/or different mix of lots of reagents. The between assay variability is ≤ 11.7%.

10. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

BIBLIOGRAPHY

1. Condemni, Jhon J. The autoimmune disease. The Journal of the American Medical Association 1987, Vo 258 n 20 2920 - 2929
2. Reichlin, M., Van Venrooij, W. J.: Autoantibodies to UNRP particles: relationship to clinical diagnosis and nephritis. Clin. exp. Immunol. 1991; 83:286-90.
3. Tan EM, et al.: The 1982 Revised Criteria for the classification of Systemic Lupus Erythematosus. Arthritis and Rheumatism 25: 1271-1277, 1982.
4. McCarty GA. Autoantibodies and their relation to rheumatic diseases. Medical Clinics of North America 70: 237-261, 1986.
5. Hardin JA. The lupus autoantigens and the pathogenesis of systemic lupus erythematosus. Arth Rheum 29: 457-460, 1986.
6. The first international standard for antibodies to double stranded DNA. Annals of the Rheumatic Disease 198; Vo 47: 740 - 746
7. Homburger HA, et al: detection of Antinuclear antibodies, comparative evaluation of enzyme immunoassay and indirect immunofluorescence methods. Arch Pathol Lab Med 122:993-999 (1998)
8. Smeenk, R. et al. Avidity of Antibodies to dsDNA. Comparison of IFT on Crithidia Luciliae, FARR assay and PEG assay. The journal of Immunology 1982 Vo 128 n.1. 73 -78

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ERROR POSSIBLE CAUSES / SUGGESTIONS

No colorimetric reaction

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation