DCM099-5  
Ed. 03/2012

## ANA Screen

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## **BIBLIOGRAFIA**

1. Antinuclear antibody The Lancet, September 15,1984:611-13
2. Reichlin,M., Van Venrooij, W. J.: Autoantibodies to UNRP particles: relationship to clinical diagnosis and nephritis.Clin.exp.Immunol. 1991; 83:286-90.
3. Tan EM, et al.: The 1982 Revised Criteria for the classification of Systemic Lupus Erythematosus. Arthritis and Rheumatism 25: 1271-1277,1982.
4. McCarty GA. Autoantibodies and their relation to rheumatic diseases. Medical Clinics of North America 70: 237-261, 1986.
5. Hardin JA. The lupus autoantigens and the pathogenesis of systemic lupus erythematosus. Arth Rheum 29: 457-460, 1986.
6. Targoff IN, Reichlin M.: Measurement of antibody to Jo-1 by ELISA and comparison to enzyme inhibitory activity. J immunol 38:2874-2882, 1987.
7. Homburger HA, et al: detection of Antinuclear antibodies, comparative evaluation of enzyme immunoassay and indirect immunofluorescence methods. Arch Pathol Lab Med 122:993-999 (1998)
8. Carey JL: Enzyme immunoassay fr antinuclear antibodies. Clin Lab Med 17:355-365 (1997)

**Ed. 03/2012**

**DCM099-5**

**DiaMetra S.r.l. Headquater:** Via Garibaldi, 18 - 20090 SEGRATE (MI) Italy

Tel. 0039-02-2139184 – 02-26921595

Fax 0039–02–2133354.

**Manufactory:** Via Pozzuolo 14, 06038 SPELLO (PG) Italy

Tel. 0039-0742–24851

Fax 0039–0742–316197

E-mail: info@diametra.com

DCM099-5  
Ed. 03/2012

## ANA Screen

for routine analysis

Semi-quantitative determination of IgG class against nuclear antigens in human serum or plasma

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO099

### INTENDED USE

ANA Screen kit is an indirect solid phase enzyme immunometric assay (ELISA) designed for semi-quantitative measurement of IgG antibodies directed against Sm (Smith), RNP/Sm, Scl-70, SS-A (Ro) (52kDa e 60kDa), SS-B (La), Jo1, U1-SmRNP, CENP-B, dsDNA, and Histones in human serum or plasma. ANA Screen kit is intended for laboratory use only.

### 1. CLINICAL SIGNIFICANCE

Autoimmune diseases are characterized by an abnormal functioning of the immune system which mistakenly attacks healthy body tissue. It results in a change of the mechanisms of cell recognition, the organism no longer distinguish the "self" from "not-self", ie the elements belonging to it and those that do not. This leads to production of antibodies that can affect individual organs or trigger systemic diseases, damaging entire functions of the individual.

The processes that lead to serious autoimmune diseases are still unknown in many ways.

Rheumatic autoimmune diseases are often associated with the presence of many autoantibodies directed against nuclear and cytoplasmic antigens. These antibodies are commonly called antinuclear antibodies (ANA) and their presence is important to reveal different rheumatic diseases such as systemic lupus erythematosus (SLE), Sjögren syndrome, scleroderma and mixed connective tissue disease (MCTD).

The presence of antinuclear antibodies (ANA) can be detected with IFA technique or with an ELISA test.

The advantages of ELISA tests are: ease of use, no interpreter influence and lower false positives sample.

### 2. PRINCIPLE

ANA Screen test is based on the binding of present antibodies with Sm, RNP/Sm, SS-A (Ro), SS-B (La), Scl-70, Jo1, U1-SmRNP, CENP-B, dsDNA and Histones antigens coated on the microplates. Any present antibodies in calibrators, controls or prediluted patient samples bind to the inner surface of the wells. After a 30 minutes incubation the microplate is washed with wash buffer for removing non-reactive serum components.

An anti-human-IgG horseradish peroxidase conjugate solution recognizes IgG class antibodies bound to the immobilized antigens. After a 30 minutes incubation any excess enzyme conjugate, which is not specifically bound is washed away with wash buffer.

A chromogenic substrate solution containing TMB is dispensed into the wells. After 15 minutes of incubation the color development is stopped by adding the stop solution. The solution turns yellow at this point. The level of colour is directly proportional to the concentration of IgG antibodies present in the original sample.

### 3. REAGENTS, MATERIALS AND INSTRUMENTATION

#### 3.1. Reagents and materials supplied in the kit

- ANA Calibrator (1 vial, 1.2 mL)  
Phosphate buffer 0.1M, NaN<sub>3</sub> < 0.1%, human serum  
CAL1 **REF DCE002/09906-0**
- Controls (2 vials, 1.2 mL each, ready to use)  
Phosphate buffer 0.1M, NaN<sub>3</sub> < 0.1%, human serum  
Negative Control **REF DCE045/09901-0**  
Positive Control **REF DCE045/09902-0**
- Sample Diluent (1 vial, 100 mL)  
Phosphate buffer 0.1M, NaN<sub>3</sub> < 0.1%  
**REF DCE053-0**
- Conjugate (1 vial, 15 mL)  
Anti h-IgG conjugated with horseradish peroxidase (HRP), BSA 0.1%, Proclin < 0.0015%  
**REF DCE002/09902-0**
- Coated Microplate  
(1 breakable microplate coated with ANA antigens)  
**REF DCE002/09903-0**
- TMB Substrate (1 vial, 15 mL)  
3,3',5,5'-tetramethylbenzidine 0.26 g/L, hydrogen peroxide 0.05%  
**REF DCE004-0**
- Stop Solution (1 vial, 15 mL)  
Sulphuric acid 0.15M **REF DCE005-0**
- 10X Conc. Wash Solution (1 vial, 50 mL)  
Phosphate buffer 0.2M, proclin < 0.0015%  
**REF DCE054-0**

#### 3.2. Reagents necessary not supplied

Distilled water.

#### 3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplate reader (450 nm)

### 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent.

Therefore, Calibrators and Controls should be handled in the same manner as potentially infectious material.

- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of Sodium Azide ( $\text{NaN}_3$ ) or Proclin 300<sup>R</sup> as preservatives. Avoid the contact with skin or mucosa.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, allow scroll through large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/ $\text{H}_2\text{O}_2$  to directed sunlight, metals or oxidants. Do not freeze the solution.

## 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinization and decontamination are effective in avoiding contamination of the conjugate; **this procedure is highly recommended when the kit is processed using analyzers which are not equipped with disposable tips.**  
For this purpose, Diametra supplies a separate decontamination reagent for cleaning needles.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate

is used, it is recommended to repeat the dose response curve in each plate

- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

## 6. PROCEDURE

### 6.1. Preparation of the Calibrator (C<sub>1</sub>)

Since no international reference preparation for ANA antibodies is available, the assay system is calibrated in relative arbitrary units. The concentration of the Calibrator is printed on the label.

### 6.2. Preparation of the Sample

For determination of ANA antibodies, human serum or plasma are the preferred sample matrixes.

**All serum and plasma samples have to be prediluted with sample diluent 1:100;** for example 10  $\mu\text{L}$  of sample may be diluted with 990  $\mu\text{L}$  of sample diluent.

The patients need not to be fasting, and no special sample preparation is necessary. Collect blood by venipuncture into vacutainers and separate serum (after clot formation) or plasma from the cells by centrifugation.

Samples may be stored refrigerated at 2-8°C for at least 5 days. For longer storage of up to six months samples should be stored frozen at -20°C. To avoid repeated thawing and freezing the samples should be aliquoted.

Neither Bilirubin nor Hemolysis have significant effect on the procedure.

The Controls are ready to use.

### 6.3. Preparation of the Wash Solution

Dilute the contents of each vial of the buffered wash solution concentrate (10X) with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C. In the concentrated wash solution is possible to observe the presence of crystals; in this case mix at room temperature until the complete dissolution of crystals; for greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

### 6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.**
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C<sub>1</sub>), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Controls	Blank
Calibrator C <sub>1</sub>	100 µL		
Controls		100 µL	
Diluted Sample		100 µL	
Incubate for 30 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells 3 times with 300 µL of diluted wash solution.			
Conjugate	100 µL	100 µL	
Incubate for 30 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells 3 times with 300 µL of diluted wash solution.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate for 15 minutes in the dark at room temperature (22-28°C).			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against Blank within 5 minutes.			

## 7. RESULTS

### 7.1. Results calculation

First of all determine the average of the absorbances for every duplicate. To obtain sample value, divide the absorbance of the sample by the absorbance of the Calibrator, then multiply for the concentration of the Calibrator, as shown below:

$$\text{Sample .Conc} = \frac{O.D.sample}{O.D.C1} \times \text{Conc .C1}$$

Reactivity is not connected in a linear way to the amount of antibodies present. Although an increase or a decrease in the concentration of antibodies results in increased or decreased responsiveness, the change is not in proportion (eg doubling the concentration of antibodies does not lead to a doubling of responsiveness). To achieve greater accuracy in the determination of antibodies it is recommended to test serial dilutions of the sample. The final dilution that is positive in the test is the antibody concentration of the patient.

## 8. REFERENCE VALUES

In a normal range study with serum samples from healthy blood donors the following ranges have been established with the anti-ANA screen test:

	ANA Screen (AU/mL)
Negativo	< 25
Positivo	≥ 25

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

Positive results should be verified concerning the entire clinical status of the patient. Also every decision for therapy should be taken on an individual patient basis. It is recommended that each laboratory establishes its own normal and pathological ranges of serum anti-ANA.

### 8.1. Specificity

Comparison test against a commercial reference kit, performed on 36 sera (28 of them positive sera and 8 negative sera) shows a 80.0 % specificity.

### 8.2. Sensitivity

Comparison test against a commercial reference kit, performed on 36 sera (28 of them positive sera and 8 negative sera) shows a 87.1 % sensitivity.

### 8.3. Detection limit

The lowest concentration of ANA that can be distinguished from the negative control is 0.69 AU/mL with a confidence limit of 95%.

### 8.4. Precision and Reproducibility

#### 8.4.1. Intra-Assay

Within run variation was determined by replicate 16 times two different sera with values in the range of calibration curve. The within assay variability is ≤ 6.6%

#### 8.4.2. Inter-Assay

Between run variation was determined by replicate measurements of two different control sera with different lots of kits and/or different mix of lots of reagents. The between assay variability is ≤ 13.3%.

## 9. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

## BIBLIOGRAPHY

1. Antinuclear antibody The Lancet, September 15,1984:611-13
2. Reichlin,M., Van Venrooij, W. J.: Autoantibodies to UNRP particles: relationship to clinical diagnosis and nephritis.Clin.exp.Immunol. 1991; 83:286-90.
3. Tan EM, et al.: The 1982 Revised Criteria for the classification of Systemic Lupus Erythematosus. Arthritis and Rheumatism 25: 1271-1277,1982.
4. McCarty GA. Autoantibodies and their relation to rheumatic diseases. Medical Clinics of North America 70: 237-261, 1986.
5. Hardin JA. The lupus autoantigens and the pathogenesis of systemic lupus erythematosus. Arth Rheum 29: 457-460, 1986.
6. Targoff IN, Reichlin M.: Measurement of antibody to Jo-1 by ELISA and comparison to enzyme inhibitory activity. J immunol 38:2874-2882, 1987.
7. Homburger HA, et al: detection of Antinuclear antibodies, comparative evaluation of enzyme immunoassay and indirect immunofluorescence methods. Arch Pathol Lab Med 122:993-999 (1998)
8. Carey JL: Enzyme immunoassay fr antinuclear antibodies. Clin Lab Med 17:355-365 (1997)

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**ERROR POSSIBLE CAUSES / SUGGESTIONS****No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

**Too low reaction (too low ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

**Too high reaction (too high ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

**Unexplainable outliers**

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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