DCM115-4  
Ed. 03/2012

# Anti Tg ELISA

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# Anti Tg ELISA

for routine analysis

Quantitative determination of antibodies against Thyroglobulin (TG) in human serum or plasma

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO115

## INTENDED USE

Indirect Immunoenzymatic colorimetric method for quantitative determination of Anti-Tg concentration in human serum or plasma.

Ani TG kit is intended for laboratory use only.

## 1. CLINICAL SIGNIFICANCE

Thyroglobulin (Tg) is a protein produced by and used entirely within the thyroid gland. Thyroglobulin is used by the thyroid gland to produce the thyroid hormones thyroxine (T4) and triiodothyronine (T3).

Thyroglobulin synthesis is stimulated at the transcriptional level by thyroid-stimulating hormone (TSH). Thyroglobulin (Tg) is a well-known target for autoantibodies occurring in thyroid autoimmunity (Graves' disease and Hashimoto's thyroiditis). Anti-Tg antibodies mostly belong to the IgG class. Low to moderate levels of anti-Tg antibodies can be found in sera of other autoimmune patients (eg systemic lupus erythematosus or Sjogren syndrome). In some cases anti-Tg positive sera may show negativity for other type of anti-thyroid antibodies - anti-TPO.

Therefore, combined determination of both types of anti-thyroid antibodies (anti-TPO + anti-Tg) provides most sensitive laboratory diagnostic tool for thyroid autoimmunity. Separately from autoimmunity, anti-Tg antibodies may develop in patients suffering from thyroid cancer. High level of anti-Tg in such patients may interfere with correct determination of serum thyroglobulin, which serves as tumour marker for therapy control in this group of patients.

## 2. PRINCIPLE

This test is based on a two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated with the antigen. Antibodies from the specimen bind to the coated antigen on the microwell surface. Unbound material is removed by the washing procedure. A second antibody directed towards human IgG and labeled with peroxidase enzyme is then added into the microwells. After the subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of a chromogen substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of specific antibodies in the specimen.

## 3. REAGENTS, MATERIALS AND INSTRUMENTATION

### 3.1. Reagents and materials supplied in the kit

- Anti-Tg Calibrators** (6 vials, 1.2 mL each)  
Phosphate buffer 0.1M, NaN<sub>3</sub> <0.1%, human serum  
CAL0 **REF DCE002/11506-0**  
CAL1 **REF DCE002/11507-0**  
CAL2 **REF DCE002/11508-0**  
CAL3 **REF DCE002/11509-0**  
CAL4 **REF DCE002/11510-0**  
CAL5 **REF DCE002/11511-0**
- Controls** (2 vials, 1.2 mL each, ready to use)  
Phosphate buffer 0.1M, NaN<sub>3</sub> <0.1%, human serum  
Negative Control **REF DCE045/11501-0**  
Positive Control **REF DCE045/11502-0**
- Sample Diluent** (1 vial, 100 mL)  
Phosphate buffer 0.1M, NaN<sub>3</sub> <0.1% **REF DCE053-0**
- Conjugate** (1 vial, 15 mL)  
Anti h-IgG conjugated with horseradish peroxidase (HRP), BSA 0.1%, Proclin <0.0015% **REF DCE002/11502-0**
- Coated Microplate**  
(1 breakable microplate coated with antigenic Thyroglobulin) **REF DCE002/11503-0**
- TMB Substrate** (1 vial, 15 mL)  
H<sub>2</sub>O<sub>2</sub>-TMB (0,26 g/L) (avoid any skin contact) **REF DCE004/-0**
- Stop Solution** (1 vial, 15 mL)  
0.15M Sulfuric acid (avoid any skin contact) **REF DCE005-0**
- 50X Conc. Wash Solution** (1 vial, 20 mL)  
NaCl 45 g/L; Tween-20 55 g/L **REF DCE006-0**

### 3.2. Reagents necessary not supplied

Distilled water.

### 3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm)

## 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.

- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators and the Controls should be handled in the same manner as potentially infectious material.
- Some reagents contain small amounts of Sodium Azide (NaN<sub>3</sub>) or Proclin 300<sup>R</sup> as preservatives. Avoid the contact with skin or mucosa.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, allow scroll through large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.

## 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinization and decontamination are effective in avoiding contamination of the conjugate; **this procedure is highly recommended when the kit is processed using analyzers which are not equipped with disposable tips.** For this purpose, Diametra supplies a separate decontamination reagent for cleaning needles.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.

- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

## 6. PROCEDURE

### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>5</sub>)

The assay system is calibrated in relative arbitrary units. The Calibrators are ready to use and have the following concentration:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
AU/mL	0	2	4	8	32	128

Once opened, the Calibrators are stable 6 months at 2-8°C.

### 6.2. Preparation of the Sample

For determination of Anti-Tg human serum or plasma are the preferred sample matrix.

**Serum samples have to be prediluted with sample diluent 1:100;** for example, 10 µL of sample may be diluted with 990 µL of sample diluent.

The patients need not to be fasting, and no special preparations are necessary. Collect blood by venipuncture into vacutainers and separate serum (after clot formation) or plasma from the cells by centrifugation.

Samples may be stored refrigerated at 2-8°C for at least 5 days. For longer storage of up to six months samples should be stored frozen at -20°C. To avoid repeated thawing and freezing the samples should be aliquoted.

Do not use microbiologically contaminated samples, as highly lipemic or hemolyzed samples.

The Controls are ready to use.

### 6.3. Preparation of the Wash Solution

Dilute the contents of each vial of the buffered wash solution concentrate (50X) with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C. In concentrated wash solution is possible to observe the presence of crystals; in this case mix at room temperature until the complete dissolution of crystals; for greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care to transfer completely the crystals, then mix until crystals are completely dissolved.

#### 6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.**
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample or Controls	Blank
Calibrator C <sub>0</sub> -C <sub>5</sub>	100 µL		
Controls		100 µL	
Diluted Sample		100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells 3 times with 300 µL of diluted wash solution			
Conjugate	100 µL	100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells 3 times with 300 µL of diluted wash solution.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate 15 minutes in the dark at room temperature (22-28°C).			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against Blank within 5 minutes.			

CAL	OD1	OD2	mean	C1	C2	mean	CV%
0	0.025	0.032	0.0285	0.045	0	0.0225	141,42
1	0.231	0.24	0.2355	1.88	2.12	2	8,49
2	0.429	0.442	0.4355	3.96	4.2	4.08	4,16
3	0.767	0.759	0.763	8.2	7.77	7.985	3,81
4	1.884	2.065	1.9745	31.9	32.36	32.13	1,01
5	3.190	3.316	3.253	128.2	127.5	127.85	0,39

#### 8. REFERENCE VALUES

In a normal range study with serum samples from healthy blood donors the following ranges have been established with the Anti-Tg tests:

	Anti Tg [AU/mL]
Normal	< 4
Elevated	≥ 4

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

Positive results should be verified concerning the entire clinical status of the patient. Also every decision for therapy should be taken individually.

It is recommended that each laboratory establishes its own normal and pathological ranges of serum Anti-Tg.

##### 8.1 Specificity

Comparison test against a commercial reference kit, performed on 66 sera (27 of them positive sera and 39 negative sera) shows a 94.9% specificity.

##### 8.2 Sensibility

Comparison test against a commercial reference kit, performed on 66 sera (27 of them positive sera and 39 negative sera) shows a 96.3% sensibility.

##### 8.3 Detection limit

The lowest concentration of anti Tg that can be distinguished from Calibrator 0 is about 0.11 AU/mL with a confidence limit of 95%.

#### 9. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

#### BIBLIOGRAPHY

- Volpè R., CRC Press , (1990)
- Volpè R., Clin. Chem. 40 2132 (1994)
- Beever K., Clin. Chem. 35 1949-54 (1989)
- Mat Clin. Chem. 40 2128 (1994)
- Chiavato L., Autoimmunity 10, 319-31 (1991)
- Degroot LJ, Thyroperoxidase Thyroid Autoimmunity 207:177 – 182 (1990)

#### 7. RESULTS

##### 7.1. Calibration curve

For Anti-Tg a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice. Smoothed-Spline Approximation and log-log coordinates are also suitable. We recommend using a Lin-Log curve.

First calculate the averaged optical densities for each calibrator well. Use lin-log graph paper and plot the averaged optical density of each calibrator versus the concentration. Draw the best fitting curve approximating the path of all calibrator points. The calibrator points may also be connected with straight line segments. The concentration of unknowns may then be estimated from the calibration curve by interpolation.

##### Typical Results (example only)

The figures below show typical results for Anti- Tg. These data are intended for illustration only and should not be used to calculate results from another run.

**Ed. 03/2012**

**DCM115-4**

**DiaMetra S.r.l. Headquarter:** Via Garibaldi, 18 – 20090

SEGRATE (MI) Italy

Tel. 0039-02-2139184 – 02-26921595

Fax 0039-02-2133354.

**Manufactory:** Via Pozzuolo 14, 06038 SPELLO (PG) Italy

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**ERROR POSSIBLE CAUSES / SUGGESTIONS****No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

**Too low reaction (too low ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

**Too high reaction (too high ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

**Unexplainable outliers**

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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