# ACE kinetic

## Angiotensin Converting Enzyme

KK-ACK	100 tests
KK-ACK2	2 x 50 tests
KK-ACKX	1200 tests

## ENGLISH

### INTENDED USE

The BÜHLMANN ACE kinetic test is intended for the direct and quantitative *in vitro* diagnostic determination of angiotensin converting enzyme (ACE) activity in serum by an enzymatic assay.

### PRINCIPLE OF THE ASSAY

ACE catalyses the conversion of angiotensin I to angiotensin II. The enzyme also mediates the cleavage of a synthetic substrate (FAPGG) into an amino acid derivative and a dipeptide. The kinetics of this cleavage reaction are measured by following the decrease in absorbance at 340 nm (1,2). The ACE kinetic method is standardized with the BÜHLMANN ACE colorimetric kit (order code: KK-ACE) according to the method described in references 3 and 4 below.

### REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity			Code	Reconstitution
	KK-ACK	KK-ACK2	KK-ACKX		
Substrate	1 vial 26 ml	2 vials 13 ml	3 vials 100 ml	B-ACK-SUB B-ACK2-SUB <sup>1</sup> B-ACKX-SUB <sup>2</sup>	Ready to use
Calibrator <sup>3</sup>	1 vial	2 vials	3 vials	B-ACK-CA	add 2 ml of sterile water
Controls <sup>4</sup> Normal and High	2 vials	2x2 vials	3x2 vials	B-ACK-CONSET	add 2 ml of sterile water

Table 1

<sup>1</sup> Order Codes for KK-ACK2.

<sup>2</sup> Order Codes for KK-ACKX.

<sup>3</sup> Lyophilized ACE Calibrator in a protein serum matrix with lot specific activity. After reconstitution leave for 15 minutes at 18-28°C and mix well before use. The Calibrator has been standardized with the BÜHLMANN ACE colorimetric kit (*cf.* above).

<sup>4</sup> Lyophilized ACE Normal and High Controls in a protein serum matrix with lot specific activity. After reconstitution leave for 15 minutes at 18-28°C and mix well before use.

### STORAGE AND SHELF LIFE OF REAGENTS

Unopened Reagents	
Stable at 2-8°C until expiration date printed on the label	
Opened / Reconstituted Reagents	
Substrate	Stable until exp. date at 2-8°C
Calibrator	Stable for 6 months at 2-8°C
ACE Controls	

Table 2

### WARNINGS AND PRECAUTIONS

All reagents of this kit contain components of human origin. Each serum donor unit used in the preparation of the kit components was tested by an FDA approved method and found negative for HBV surface antigen, so as for HCV and HIV1/2 antibodies. Although these methods are highly accurate, there is no guarantee that this material cannot transmit Hepatitis or AIDS. *Therefore, all patient specimens and kit components should be handled as if capable of transmitting infections.* All products containing human source material should be handled in accordance with good laboratory practice using appropriate precautions.

Calibrator and ACE Controls contain sodium azide as a preservative in a concentration less than 0.01 mg/ml

### MATERIALS REQUIRED BUT NOT PROVIDED

- Precision pipettes: 25, 100, 150, 250 µl and 2 ml
- Plastic test tubes for sample or calibrator incubation
- Vortex mixer
- Waterbath for incubation set at 37°C

- Spectrophotometer with temperature controlled cuvette holders for incubation at 37°C and for measurement of absorbance at 340 nm
- Clinical chemistry analyzer with 340 nm filter (optional)

### SPECIMEN COLLECTION AND STORAGE

Since EDTA inhibits ACE activity, serum specimen should be used for the determination of ACE activity:

- Collect sufficient blood (at least 0.5 ml) by venipuncture into an appropriate tube without anticoagulant.
- Let the blood stand for 2-3 hours at room temperature, then centrifuge at 4°C and 1000 x g; collect the serum.
- Freeze the serum specimen at -20°C if not assayed within 5 days. ACE activity in sterile serum is stable for up to 30 days at 2-8°C and 6 months at -20°C.

To avoid lipemic sera, the blood sample should be taken from fasting patient. Due to interference with the photometric determination lipemic sera must be pretreated either with an Ultracentrifuge or with LipoClear from StatSpin Inc. ([www.statspin.com](http://www.statspin.com)). Icteric or hemolytic sera can not be used for the ACE activity determination.

### ASSAY PROCEDURE

**Assay Protocol** (For automatic procedure *cf.* next chapter below)

Allow the substrate to come to room temperature. Avoid heating of the substrate. Pretreat all lipemic sera according to the procedure as described above.

1. Label two plastic tubes for each calibrator and control serum. Label additional tubes for each patient sample.
2. Pipet 25 µl of Calibrator, Control serum and patient Sample, respectively, into the corresponding tubes.

NOTE: For each sample follow individually the procedure outlined below (according to the capacity of the photometer used):

3. Add 250 µl of substrate to each sample and vortex thoroughly.
4. Incubate for 5 minutes at 37°C in a water bath.
5. Set your photometer to zero absorbance with distilled water.
6. Transfer sample into a microcuvette and measure absorbance at 37°C and 340 nm twice in a time interval of exactly 10 minutes.

### RESULTS AND AUTOMATION

#### Standardization

The values of this assay are identical to the values obtained with the standard BÜHLMANN ACE colorimetric assay kit (order code: KK-ACE). This standardization has the advantage that results obtained with either the colorimetric method or the radioenzymatic method are intercomparable.

#### Automation with Chemistry Analyzers

ACE kinetic can be performed on every open chemistry analyzer. For instructions on programming and operation please refer to the "Operators Manual" of the Instrument.

NOTE: Each laboratory should establish its own instrument settings for the individual requirements.

PARAMETER LISTINGS FOR SEVERAL CHEMICAL ANALYZERS ARE AVAILABLE UPON REQUEST.

#### CALCULATION

Calculate the corresponding enzyme activity ( $E_x$ ) of each unknown sample by dividing the difference in absorbance of the individual sample ( $\Delta A_x$ ) through the mean of the absorbance difference of the calibrator vial ( $\Delta A_c$ ) and multiply by the enzyme activity ( $E_c$ ) indicated on the data sheet (for an example see Table 11):

$$E_x = \frac{\Delta A_x}{\Delta A_c} \times E_c$$

**Definition:** One unit of ACE activity is defined as the amount of enzyme required to release one  $\mu\text{mol}$  of Hippuric Acid per minute and per liter of serum at 37°C as determined by the colorimetric assay:

$$1 \text{ ACE unit} = \frac{1 \mu\text{mol hippuric acid}}{\text{min} \times L}$$

**QUALITY CONTROL**

The values of the Normal and High Controls provided with the kit must be within the lot specific range indicated on the corresponding data sheet. If not, the assay has to be repeated.

It is good laboratory practice to record the following data for each assay: kit lot number, reconstitution dates of kit components used, concentration value of calibrator and of controls, concentration values of internal serum pool.

**PERFORMANCE CHARACTERISTICS**

**Intra-Assay Precision: 2.7%.** The intra-assay precision has been determined by measuring three different serum samples 20 times in a single test run (cf. Table 12).

**Inter-Assay Precision: 8.1%.** The inter-assay precision has been determined by measuring three control serum samples in 20 different runs (cf. Table 13).

**Sensitivity:** The minimum detectable dose of ACE was determined to be **12 ACE units**.

**Dilution Linearity: r = 0.998.** Two elevated serum samples have been diluted with physiological NaCl solution (cf. Table 14).

**Spiking Recovery: 99.8%.** Two different serum samples have been spiked with increasing amounts of ACE and analyzed according to the assay procedure (cf. Table 15).

**Specificity:** Inhibition by its natural substrate Ang I and by EDTA and H-Val-Trp-OH. A serum sample of defined ACE activity has been inhibited and measured again (cf. Table 16).

**Linearity:** The kinetic method is linear from **12 to 150 ACE units**.

**EXPECTED VALUES**

From a group of 93 adult (age 20-70 years) blood donors from Switzerland, the normal serum range (Mean + 2SD) has been determined to be:

**12 – 68 ACE units**

Serum samples from 50 sarcoidosis patients were within a range of 45-135 units (5)

Serum ACE levels in individuals under 18 years of age are substantially higher and more variable than in adults. Those should be established separately (6).

**METHOD COMPARISON**

The comparison of 66 sera from apparently healthy donors (age: 20 – 70 years) were analyzed using the presented ACE kinetic assay as well as the ACE colorimetric and the ACE direct radioenzymatic (REA) kits from BÜHLMANN LABORATORIES AG (order codes: KK-ACE and RK-ACD). The subsequent linear regression analysis yielded the following statistics:

ACE kinetic = 1.098 \* ACE colorimetric – 7.83 ACE units;  
r = 0.964; R<sup>2</sup> = 0.868

ACE kinetic = 1.449 \* ACE direct – 17.22 ACE units;  
r = 0.938; R<sup>2</sup> = 0.880

Means: 42.1 ACE units (ACE kinetic)  
45.5 ACE units (ACE colorimetric)  
40.9 ACE units (ACE direct)

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