

# OnSite™ HSV-1 IgG/IgM Rapid Test

**REF R0203C** 

### INTENDED USE

The OnSite HSV-1 IgG/IgM Rapid Test is a lateral flow chromatographic immunoassay for the qualitative detection and differentiation of antibodies (IgG and IgM) to herpes simplex virus 1 (HSV-1) in human serum, plasma or whole blood. It is intended to be used by professionals as a screening test and as an aid in the diagnosis of infection with HSV-1. Any reactive result with the OnSite HSV-1 IgG/IgM Rapid test must be confirmed with alternative testing method(s) and clinical findings.

### SUMMARY AND EXPLANATION OF THE TEST

Herpes simplex viruses are two types of DNA viruses of the Herpesviridae family, herpes simplex virus-1 (HSV-1) and herpes simplex virus-2 (HSV-2). HSV-1 is generally acquired during childhood via non-sexual contact and affects mainly the orofacial area. HSV-2 is nearly always sexually transmitted and is the main cause of genital herpes. HSV-1 and HSV-2 can infect both genital and orofacial areas<sup>1</sup>. Up to 50% of first-episode cases of genital herpes are caused by HSV-1, but recurrences are much less frequent for genital HSV-1 infection than genital HSV-2 infection<sup>2</sup>. HSV subclinical viral shedding are less frequent for genital HSV-1 than genital HSV-2<sup>2</sup>. Genital HSV infection has also been associated with increased risk for sexual transmission of HIV<sup>3</sup>. After primary infection, these viruses persist in a latent state for life<sup>1</sup>.

One of the biggest risks associated with HSV is neonatal transmission<sup>1</sup>. The rate of neonatal transmission is higher in mothers with genital HSV-1 than in those with genital HSV-2<sup>4</sup>. Eighty-five to ninety percent of neonatal transmission occurs at the time of delivery with only 5% of infections occurring intrauterine<sup>5</sup>. Clinical manifestations of neonatal infection with HSV range from local lesions of the skin, mouth, eye or central nervous system to severe, widespread dissemination involving visceral organs and potentially death<sup>1</sup>.


Serology is an effective means of diagnosing HSV because the manifestation of symptoms is transient and the infection is often undiagnosed<sup>1</sup>. IgM anti-HSV can be detected 9-10 days after exposure and last for 7-14 days, although it may remain detectable for up to 6 weeks<sup>6</sup>. IgM anti-HSV is often associated with primary infection but may be detectable during recurrence of the disease<sup>6</sup>. IgG anti-HSV can be detected 21-28 days post exposure and detectable titers typically remain for life<sup>6</sup>. Detection of IgM anti-HSV in the absence of IgG anti-HSV can be an effective tool in detecting early stages of HSV infection and as an indicator of potential primary infection.

HSV-1 and HSV-2 infections have different prognoses. Type-specific serological diagnosis is beneficial and can be achieved by using glycoprotein G1 and glycoprotein G2 as recommended by the CDC<sup>7</sup>.

The OnSite HSV-1 IgG/IgM Rapid Test uses HSV-1 glycoprotein G1 for the specific detection and differentiation of IgG and IgM antibodies to HSV-1 in serum, plasma and whole blood. The test can be performed in 10 minutes by minimally skilled personnel without the use of laboratory equipment.

### TEST PRINCIPLE

The OnSite HSV-1 IgG/IgM Rapid Test is a lateral flow chromatographic immunoassay. The test cassette consists of: 1) a burgundy colored conjugate pad containing HSV-1 type specific glycoprotein G1 antigens conjugated with colloidal gold (HSV-1 conjugates) and a control antibody conjugated with colloidal gold, 2) a nitrocellulose membrane strip containing two test lines (G and M lines) and a control line (C line). The G line is pre-coated with mouse anti-human IgG for detection of IgG anti-HSV-1, the M line is pre-coated with mouse anti-human IgM for detection of IgM anti-HSV-1, and the C line is pre-coated with a control line antibody.



When an adequate volume of test specimen and sample diluent are dispensed into the sample well, the specimen migrates by capillary action across the cassette. IgG anti-HSV-1, if present in the specimen, will bind to the HSV-1 conjugates. The immunocomplex is then captured on the membrane by the pre-coated mouse anti-human IgG forming a burgundy colored G line, indicating an HSV-1 IgG positive test result. IgM anti-HSV-1, if present in the specimen, will bind to the HSV-1 conjugates. The immunocomplex is then captured on the membrane by the pre-coated mouse anti-human IgM forming a burgundy colored M line, indicating an HSV-1 IgM positive test result.

Absence of any test lines (G or M) suggests a negative result. The test contains an internal control (C line) which should exhibit a burgundy colored line of the immunocomplex of the control antibodies, regardless of color development on the test lines (G and M). If no control line (C line) develops, the test result is invalid and the specimen must be retested with another device.

### REAGENTS AND MATERIALS PROVIDED

- Individually sealed foil pouch containing:
  - One cassette device
  - One desiccant
- 10 µL capillary tubes
- Sample diluent (REF SB-R0203, 5 mL/bottle)
- One package insert (instruction for use)

### MATERIALS MAY BE REQUIRED AND NOT PROVIDED

- Positive control
- Negative control

### MATERIALS REQUIRED BUT NOT PROVIDED

- Clock or Timer
- Lancing device for whole blood test

### WARNINGS AND PRECAUTIONS

#### For in vitro Diagnostic Use

- This package insert must be read completely before performing the test. Failure to follow the insert may lead to inaccurate test results.
- Do not open the sealed pouch until ready to conduct the assay.

- Do not use expired devices or components.
- Bring all reagents to room temperature (15-30°C) before use.
- Do not use components from another test kit to substitute for components of this kit.
- Do not use hemolyzed blood specimens for testing.
- Wear protective clothing and disposable gloves while handling the kit reagents and clinical specimens. Wash hands thoroughly after performing the test.
- Users of this test should follow the US CDC Universal Precautions for prevention of transmission of HIV, HBV and other blood-borne pathogens.
- Do not smoke, drink or eat in areas where specimens or kit reagents are being handled.
- Dispose of all specimens and materials used to perform the test as bio-hazardous waste.
- Handle negative and positive controls in the same manner as patient specimens.
- The test result should be read 10 minutes after a specimen is applied to the sample well of the device. Reading the test result after 15 minutes may give erroneous results.
- Do not perform the test in a room with strong air flow, e.g. an electric fan or strong air conditioning.

### REAGENT PREPARATION AND STORAGE INSTRUCTIONS

All reagents are ready to use as supplied. Store unused test devices unopened at 2-30°C. If stored at 2-8°C, ensure that the test device is brought to room temperature before opening. The test device is stable through the expiration date printed on the sealed pouch. Do not freeze the kit or expose the kit to temperatures above 30°C.

### SPECIMEN COLLECTION AND HANDLING

Consider any materials of human origin as infectious and handle them using standard bio-safety procedures.

#### Plasma

- Collect blood specimen into a lavender, blue or green top collection tube (containing EDTA, citrate or heparin, respectively, in Vacutainer®) by venipuncture.
- Separate the plasma by centrifugation.
- Carefully withdraw the plasma into new pre-labeled tube.

#### Serum

- Collect blood specimen into a red top collection tube (containing no anticoagulants in Vacutainer®) by venipuncture.
- Allow the blood to clot.
- Separate the serum by centrifugation.
- Carefully withdraw the serum into a new pre-labeled tube.

Test specimens as soon as possible after collecting. Store specimens at 2-8°C, if not tested immediately, for up to 5 days. The specimens should be frozen at -20°C for longer storage.

Avoid multiple freeze-thaw cycles. Prior to testing, bring frozen specimens to room temperature slowly and mix gently. Specimens containing visible particulate matter should be clarified by centrifugation before testing.

Do not use specimens demonstrating gross lipemia, gross hemolysis or turbidity in order to avoid interference with result interpretation.



#### Whole Blood

Drops of whole blood can be obtained by either finger tip puncture or venipuncture. Collect blood specimen into a lavender, blue or green top collection tube (containing EDTA, citrate or heparin, respectively in Vacutainer®). Do not use hemolyzed blood for testing.

Whole blood specimens should be stored in refrigeration (2-8°C) if not tested immediately. The specimens must be tested within 24 hours of collection.

### ASSAY PROCEDURE

- Bring the specimen and test components to room temperature if refrigerated or frozen. Once the specimen is thawed, mix well prior to performing the assay.
- When ready to test, open the pouch at the notch and remove the device. Place the test device on a clean, flat surface.
- Be sure to label the device with the specimen's ID number.
- Fill the capillary tube with specimen not exceeding the specimen line as shown in the images below. The volume of specimen is approximately 10 µL.  
**For better precision, transfer specimen using a pipette capable of delivering a 10 µL volume.**



- Set up the timer.
- Result should be read in 10 minutes. Positive results may be visible as soon as 1 minute.


**Do not read the result after 15 minutes. To avoid confusion, discard the test device after interpreting the result.**

**QUALITY CONTROL**

- Internal Control**  
This test contains a built-in control feature, the C line. The C line develops after adding the specimen and the sample diluent. If the C line does not develop, review the entire procedure and repeat the test with a new device.
- External Control**  
Good Laboratory Practice recommends using external controls, positive and negative, to assure the proper performance of the assay, particularly under the following circumstances:
  - New operator uses the kit, prior to performing the testing of the specimens.
  - A new lot of test kits is used.
  - A new shipment of test kits is used.
  - The temperature used during storage of the kits falls outside of 2-30°C.
  - The temperature of the test area falls outside of 15-30°C.
  - To verify a higher than expected frequency of positive or negative results.
  - To investigate the cause of repeated invalid results.

**INTERPRETATION OF ASSAY RESULT**

- NEGATIVE RESULT**  
If only the C line develops, the test indicates that anti-HSV-1 antibodies are not detected in the specimen. The result is negative or non-reactive.



- POSITIVE RESULT**
  - In addition to the presence of the C line, if only the G line develops, the test result indicates the presence of IgG anti-HSV-1; the result is HSV-1 IgG positive or reactive.



- In addition to the presence of the C line, if only the M line develops, the test indicates the presence of IgM anti-HSV-1. The result is HSV-1 IgM positive or reactive.



- In addition to the presence of C line, if both the G and M lines develop, the test indicates the presence of IgG anti-HSV-1 and IgM anti-HSV-1. The result is HSV-1 IgG and HSV-1 IgM positive or reactive.



*Samples with positive results should be confirmed with alternative testing method(s) and clinical findings before a diagnosis is made.*

- INVALID**  
If no C line develops, the assay is invalid regardless of any color development on the test lines (G and M) as indicated below. Repeat the assay with a new device.



**PERFORMANCE CHARACTERISTICS**

- Accuracy of IgG Detection**  
A total of 227 specimens were collected and tested with the OnSite HSV-1 IgG/IgM Rapid Test and by a commercial IgG anti-HSV-1 ELISA. Comparison for all subjects is shown in the following table:

Reference	OnSite HSV-1 IgG/IgM Rapid Test		Total
	Positive	Negative	
Positive	174	18	192
Negative	3	32	35
<b>Total</b>	<b>177</b>	<b>50</b>	<b>227</b>

Relative Sensitivity: 90.6%, Relative Specificity: 91.4%, Overall Agreement: 90.7%

- Accuracy of IgM Detection**  
A total of 107 specimens were collected and tested with the OnSite HSV-1 IgG/IgM Rapid Test and by a commercial IgM ELISA. The overall agreement between the two products was 85.0%.

- Positive Rate on Random Clinical Specimens**  
Ten thousand random, clinical specimens were tested with the OnSite HSV-1 IgG/IgM Rapid Test. The positive rate was 92.8% for IgG anti-HSV-1 and 4.9% for IgM anti-HSV-1.

- Cross-Reactivity**  
No false positive IgG and IgM anti-HSV-1 results were observed on 3-10 specimens from the following disease states or special conditions, respectively:
 

<i>T. pallidum</i>	<i>H. pylori</i>	Dengue	Malaria	Typhoid
Toxoplasma	Rubella	CMV	HSV-2	hCG
HAV	HBV	HCV	HEV	HIV
TB	ANA	HAMA	RF (up to 1,000 IU/mL)	

- Interference**  
Common substances (such as pain and fever medication and blood components) may affect the performance of the OnSite HSV-1 IgG/IgM Rapid Test. This was studied by spiking these substances into negative, IgG positive and IgM positive specimens, respectively. The results demonstrate that at the concentrations tested, the substances studied do not affect the performance of the OnSite HSV-1 IgG/IgM Rapid Test.

List of potentially interfering substances and concentrations tested:

1. Albumin	60 g/L	6. Hemoglobin	2 g/L
2. Bilirubin	20 mg/dL	7. Heparin	3,000 U/L
3. Creatinine	442 µmol/L	8. Salicylic acid	4.24 mmol/L
4. EDTA	3.4 µmol/L	9. Sodium citrate	3.8%
5. Glucose	55 mmol/L		

**EXPECTED VALUES**

In non-high-risk populations, HSV-1 prevalence tends to increase with age, with acquisition occurring primarily in childhood and adolescence. Prevalence commonly reaches 40% by the age of 15, before increasing to 60-90% in older adults. In a given population and age group, HSV-1 prevalence is nearly always higher than HSV-2 prevalence<sup>8</sup>. Although genital herpes is primarily associated with HSV-2, an increasing proportion of genital herpes is caused by HSV-1, particularly in Europe<sup>9</sup>. Clinical studies report detection of HSV-1 IgM and IgG in 5.9% and 93.2% of patients, respectively<sup>10</sup>.

**LIMITATIONS OF TEST**

- The Assay Procedure and the Interpretation of Assay Result sections must be followed closely when testing for the presence of antibodies to HSV-1 in serum, plasma or whole blood from individual subjects. Failure to follow the procedure may lead to inaccurate results.
- The OnSite HSV-1 IgG/IgM Rapid Test is limited to the qualitative detection of antibodies to HSV-1 in human serum, plasma or whole blood. The intensity of the test line does not have linear correlation with the titer of HSV-1 antibody in the specimen.
- A negative or non-reactive test result does not preclude the possibility of exposure to or infection with HSV-1. A negative or non-reactive result can occur if the titer of HSV-1 antibody present in the specimen is below the level detectable by the assay or if HSV-1 antibody was not present during the stage of disease in which the sample was collected.
- A negative result does not rule out an infection with HSV-1. Samples collected too early in the course of an infection may not have detectable levels of IgM.
- Infection may progress rapidly. If the symptom persists, while the result from OnSite HSV-1 IgG/IgM Rapid Test is negative or non-reactive, it is recommended to re-test the patient a few days later or test with an alternative test method.
- The OnSite HSV-1 IgG/IgM Rapid Test has not been validated on specimens from neonates.
- Specimens from patients with infectious mononucleosis or high titers of heterophile antibodies, rheumatoid factor (>1,000 IU/mL) may affect expected results.
- Results obtained with the OnSite HSV-1 IgG/IgM Rapid Test should only be interpreted in conjunction with other diagnostic procedures and clinical findings.

**REFERENCES**

- Kimberlin DW. Neonatal herpes simplex infection. Clin Microbiol Rev 2004. 17(1):1-13.
- LeGoff J, Péré H, Bélec L. Diagnosis of genital herpes simplex virus infection in the clinical laboratory. Virol J. 2014. 11(1):83.
- Laderman EI, Whitworth E, Dumauld E, et al. Rapid, sensitive, and specific lateral-flow immunochromatographic point-of-care device for detection of herpes simplex virus type 2-specific immunoglobulin G antibodies in serum and whole blood. Clin Vaccine Immunol 2008. 15(1):159-63.
- Brown ZA, Wald A, Morrow RA, et al. Effect of serologic status and cesarean delivery on transmission rates of herpes simplex virus from mother to infant. J Am Med Assoc 2003. 289(2):203-9.
- Starface G, Selmin A, Zanardo V, et al. Herpes simplex virus infection in pregnancy. Infect Dis Obstet Gynecol 2012.
- Page J, Taylor J, Tideman RL, et al. Is HSV serology useful for the management of first episode genital herpes? Sex Transm Infect 2003. 79(4):276-9.
- Workowski KA, Levine WC. Sexually Transmitted Diseases Treatment Guidelines. CDC 2002. 51(RR06):1-80.
- Smith JS, Robinson NJ: Age specific prevalence of infection with herpes simplex virus type 2 and 1: a global review. J Infection Dis 2002. 86(Suppl 1):S3-82.
- Azwa A, Barton SE. Aspects of herpes simplex virus: a clinical review. J Fam Plan Reprod H 2009. 35(4):237-42.
- Obeid OE. Prevalence of herpes simplex virus types 1 and 2 and associated sociodemographic variables in pregnant women attending King Fahd hospital of the university. J Family Community Med 2007. 14(1):3-7.

**Index of CE Symbols**

	Consult instructions for use		For <i>in vitro</i> diagnostic use only		Use by
	Catalog #		Lot Number		Tests per kit
	Store between 2-30°C		Authorized Representative		Do not reuse
	Manufacturer		Date of manufacture		

**CTK Biotech, Inc.**  
10110 Mesa Rim Road  
San Diego, CA 92121, USA  
Tel: 858-457-8698  
Fax: 858-535-1739  
E-mail: info@ctkbiotech.com

**EC REP** **MDSS GmbH**  
Schiffgraben 41, 30175 Hannover, Germany  
PI-R0203C Rev. A  
Effective date: 2015-07-24  
English Version

For Export Only. Not For Resale In The USA.