Instructions for use
Glutamate ELISA

Please use only the valid version of the Instructions for Use provided with the kit

1. Introduction

1.1 Intended use and principle of the test

Enzyme Immunoassay for the quantitative determination of L-glutamate in urine, to evaluate L-glutamate homeostasis.

The determination of L-glutamate in urine is helpful for the determination of neurostress.

After extraction and derivatisation Glutamate is quantitatively determined by ELISA.

The competitive ELISA uses the microtiter plate format. The antigen is bound to the solid phase of the microtiter plate. The derivatized analyte concentrations in the standards, controls and samples and the solid phase bound analyte compete for a fixed number of antibody binding sites. When the system is in equilibrium, free antigen and free antigen-antibody complexes are removed by washing. The antibody bound to the solid phase is detected by an anti-rabbit IgG-peroxidase conjugate using TMB as a substrate. The reaction is monitored at 450 nm.

Quantification of unknown samples is achieved by comparing their absorbance with a standard curve prepared with known standards.

1.2 Clinical application

Glutamate, also known as glutamic acid, is one of the most important excitatory neurotransmitter in the central nervous system (CNS). It is released presynaptically and it binds postsynaptically to specific receptors for glutamate. The enzyme glutamic acid decarboxylase is able to convert L-glutamate in the CNS by decarboxylation to γ -aminobutyric acid (GABA), which acts as an inhibitory neurotransmitter.

Many publications postulate, that the determination of L-glutamate in urine is helpful for the determination of neurostress. The collective term „neurostress“ refers to many physical and psychological complaints caused due to our modern way of life, unfavourable environmental conditions, poor diet, medications, occupational and social stress, sleep deprivation, overstimulation or genetic predisposition. The resulting symptoms are burnout, depression, insomnia, chronic fatigue syndrome, fibromyalgia, multiple chemical sensitivity and other chronic pathologies. Furthermore, many publications describe an impaired or dysregulated glutamic acid balance in relation to previously listed symptoms. Based on this many laboratories offer the determination of specific neurostress profiles including, among others, the determination of glutamate. In this case, glutamate is detected primarily in the second morning urine by several methods. To generate a profound neurostress profile, also other analytes (serotonin, norepinephrine, dopamine, epinephrine, melatonin, DHEA and cortisol) should be determined.

2. Procedural cautions, guidelines, warnings and limitations

2.1 Procedural cautions, guidelines and warnings

- (1) This kit is intended for professional use only. Users should have a thorough understanding of this protocol for the successful use of this kit. Only the test instruction provided with the kit is valid and has to be used to run the assay. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
- (2) This assay was validated for a certain type of sample as indicated in Intended Use (please refer to Chapter 1). Any off-label use of this kit is in the responsibility of the user and the manufacturer cannot be held liable.
- (3) The principles of Good Laboratory Practice (GLP) have to be followed.
- (4) In order to reduce exposure to potentially harmful substances, wear lab coats, disposable protective gloves and protective glasses where necessary.
- (5) All kit reagents and specimens should be brought to room temperature and mixed gently but thoroughly before use. Avoid repeated freezing and thawing of reagents and specimens.
- (6) For dilution or reconstitution purposes, use deionized, distilled, or ultra-pure water.
- (7) The microplate contains snap-off strips. Unused wells must be stored at 2 – 8 °C in the sealed foil pouch with desiccant and used in the frame provided.
- (8) Duplicate determination of sample is highly recommended to be able to identify potential pipetting errors.
- (9) Once the test has been started, all steps should be completed without interruption. Make sure that the required reagents, materials and devices are prepared ready at the appropriate time.
- (10) Incubation times do influence the results. All wells should be handled in the same order and time intervals.
- (11) To avoid cross-contamination of reagents, use new disposable pipette tips for dispensing each reagent, sample, standard and control.
- (12) A standard curve must be established for each run.
- (13) The controls should be included in each run and fall within established confidence limits. The confidence limits are listed in the QC-Report provided with the kit.
- (14) Do not mix kit components with different lot numbers within a test and do not use reagents beyond expiry date as shown on the kit labels.
- (15) Avoid contact with Stop Solution containing 0.25 M H_2SO_4 . It may cause skin irritation and burns. In case of contact with eyes or skin, rinse off immediately with water.

- (16) TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them.
- (17) For information on hazardous substances included in the kit please refer to Safety Data Sheet (SDS). The Safety Data Sheet for this product is made available directly on the website of the manufacturer or upon request.
- (18) Kit reagents must be regarded as hazardous waste and disposed of according to national regulations.
- (19) The expected reference values reported in this test instruction are only indicative. It is recommended that each laboratory establishes its own reference intervals.
- (20) In case of any severe damage to the test kit or components, the manufacturer has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components must not be used for a test run. They must be stored properly until the manufacturer decides what to do with them. If it is decided that they are no longer suitable for measurements, they must be disposed of in accordance with national regulations.

2.2 Limitations

Any inappropriate handling of samples or modification of this test might influence the results.

2.2.1 Interfering substances

Urine

Please note the sample preparation stabilization of the urine sample! It cannot be excluded that high acid concentrations lead to incorrect results. Up to 20 µl 6 N HCl per 1 ml urine no influence on the results was observed.

2.2.2 Drug interferences

There are no known substances (drugs, food) which ingestion interferes with the measurement of glutamate level in the sample.

2.2.3 High-Dose-Hook effect

No hook effect was observed in this test.

3. Storage and stability

Store the unopened reagents at 2 – 8 °C until expiration date. Do not use components beyond the expiry date indicated on the kit labels. Once opened the reagents are stable for 2 months when stored at 2 – 8 °C. Once the resealable pouch has been opened, care should be taken to close it tightly with desiccant again.

4. Materials

4.1 Contents of the kit

BA D-0090 **FOILS** **Adhesive Foil** – Ready to use

Contents: Adhesive Foils in a resealable pouch

Volume: 1 x 4 foils

BA D-0024 **REAC-PLATE** **Reaction Plate** – Ready to use

Contents: 1 x 96 well plate, empty in a resealable pouch

BA E-2442 **EXTRACT-PLATE 48** **Extraction Plate** – Ready to use

Contents: 2 x 48 well plate, precoated with cation exchanger in a resealable pouch

BA E-0030 **WASH-CONC 50x** **Wash Buffer Concentrate** – Concentrated 50x

Contents: Buffer with a non-ionic detergent and physiological pH

Volume: 1 x 20 ml/vial, light purple cap

BA E-0040 **CONJUGATE** **Enzyme Conjugate** – Ready to use

Contents: Goat anti-rabbit immunoglobulins conjugated with peroxidase

Volume: 1 x 12 ml/vial, red cap

BA E-0055 **SUBSTRATE** **Substrate** – Ready to use

Contents: Chromogenic substrate containing tetramethylbenzidine, substrate buffer and hydrogen peroxide

Volume: 1 x 12 ml/black vial, black cap

BA E-0080 **STOP-SOLN** **Stop Solution** – Ready to use

Contents: 0.25 M sulfuric acid
 Volume: 1 x 12 ml/vial, light grey cap
 Hazards identification:

H290 May be corrosive to metals.

BA E-2431 **GLUT** **Glutamate Microtiter Strips** – Ready to use

Contents: 1 x 96 well (12x8) antigen precoated microwell plate in a resealable foil pouch with desiccant

BA E-2410 **AS GLUT** **Glutamate Antiserum** – Ready to use

Contents: Rabbit anti-glutamate antibody, blue coloured
 Volume: 1 x 6 ml/vial, blue cap

BA E-2413 **ASSAY-BUFF** **Assay Buffer** – Ready to use

Contents: Buffer with alkaline pH
 Volume: 1 x 20 ml/vial, yellow cap

BA E-2428 **EQUA-REAG** **Equalizing Reagent** – Lyophilized

Contents: Lyophilized protein
 Volume: 1 vial, brown cap

Standards and Controls – Ready to use

Cat. no.	Component	Colour/Cap	Concentration µg/ml	Concentration µmol/l	Volume/ Vial
BA E-2401	STANDARD A	white	0	0	4 ml
BA E-2402	STANDARD B	light yellow	0.6	4.08	4 ml
BA E-2403	STANDARD C	orange	2	13.6	4 ml
BA E-2404	STANDARD D	dark blue	6	40.8	4 ml
BA E-2405	STANDARD E	light grey	20	136	4 ml
BA E-2406	STANDARD F	black	60	408	4 ml
BA E-2451	CONTROL 1	light green	Refer to QC-Report for expected value and acceptable range!		4 ml
BA E-2452	CONTROL 2	dark red			4 ml

Conversion: Glutamate (µg/ml) x 6.8 = Glutamate (µmol/l)

Contents: Acidic buffer with non-mercury preservative, spiked with defined quantity of Glutamate

BA E-2446 **D-REAGENT** **D-Reagent** – Ready to use

Contents: Crosslinking agent in dimethylsulfoxide
 Volume: 1 x 3 ml/vial, white cap

Hazards identification:

H317 May cause an allergic skin reaction.

BA E-2458 **Q-BUFFER** **Q-Buffer** – Ready to use

Contents: TRIS buffer
 Volume: 1 x 20 ml/vial, white cap

BA E-2460 **DILUENT** **Diluent** – Ready to use

Contents: Buffer with sodium acetate
 Volume: 1 x 20 ml/vial, dark green cap

BA E-2787**NaOH****NaOH** – Ready to use

Contents: Sodium hydroxide solution

Volume: 1 x 2 ml/vial, purple cap

Hazards
identification:

H290 May be corrosive to metals.

H315 Causes skin irritation.

H319 Causes serious eye irritation.

4.2 Additional materials and equipment required but not provided in the kit

- Calibrated precision pipettes to dispense volumes between 10 – 100 µl; 12.5 ml
- Microtiter plate shaker (shaking amplitude 3 mm; approx. 600 rpm)
- ELISA reader capable of reading absorbance at 450 nm and if possible 620 – 650 nm
- Shaker (shaking amplitude 3 mm; approx. 600 rpm)
- Absorbent material (paper towel)
- Vortex mixer
- Water (deionized, distilled, or ultra-pure)

5. Sample collection and storage**Urine**

Spontaneous urine (second morning urine) stabilized with 10 µl 6 M HCl per 1 ml of urine sample can be used.

Storage: up to 6 hours (18 – 25 °C); up to 14 days (2 – 8 °C); up to 6 months (< -15 °C).

Repeated freezing and thawing should be avoided. Avoid exposure to direct sunlight.

6. Test procedure

Allow all reagents and samples to reach room temperature and mix thoroughly by gentle inversion before use. Duplicate determinations are recommended. It is recommended to number the strips of the microwell plate before usage to avoid any mix-up.

The binding of the antisera and of the enzyme conjugate and the activity of the enzyme are temperature dependent. The higher the temperature, the higher the absorption values will be. Varying incubation times will have similar influences on the absorbance. The optimal temperature during the Enzyme Immunoassay is between 20 – 25 °C.

During the overnight incubation at 2 – 8 °C with the antiserum, the temperature should be uniform all over the ELISA plate to avoid any drift and edge-effect.

6.1 Preparation of reagents**Wash Buffer**

Dilute the 20 ml Wash Buffer Concentrate with water (deionized, distilled, or ultra-pure) to a final volume of 1000 ml.

Storage: 2 months at 2 – 8 °C

Equalizing ReagentReconstitute the Equalizing Reagent with **12.5 ml** of **Assay Buffer**.

Reconstituted Equalizing Reagent which is not used immediately has to be stored in aliquots for max 2 months at -20 °C and may be thawed only once.

D-Reagent

The D-Reagent has a freezing point of 18.5 °C. To ensure that the D-Reagent is liquid when being used, it must be ensured that the D-Reagent has reached room temperature and forms a homogeneous, crystal-free solution.

Glutamate Microtiter Strips

In rare cases residues of the blocking and stabilizing reagent can be seen in the wells as small, white dots or lines. These residues do not influence the quality of the product.

Extraction Plate

In rare cases residues of the cation exchanger can be seen in the wells as small, black dots or lines. These residues do not influence the quality of the product.

6.2 Preparation of samples

Extraction

1.	Pipette 100 µl of the standards, controls and urine samples into the appropriate wells of the Extraction Plate .
2.	Add 100 µl of the Diluent to all wells. Cover plate with Adhesive Foil and shake for 10 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm).
3.	Use 25 µl for the subsequent derivatization !

Derivatization

1.	Pipette 25 µl of the extracted standards, controls and urine samples into the appropriate wells of the Reaction Plate .
2.	Pipette 10 µl of NaOH into all wells.
3.	Pipette 50 µl of the Equalizing Reagent into all wells.
4.	Pipette 10 µl of the D-Reagent into all wells.
5.	Cover plate with Adhesive Foil and shake for 2 h at RT (20 – 25 °C) on a shaker (approx. 600 rpm).
6.	Pipette 75 µl of the Q-Buffer into all wells.
7.	Shake for 10 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm).
8.	Use 25 µl for the ELISA !

6.3 Glutamate ELISA

1.	Pipette 25 µl of the prepared standards, controls and urine samples into the appropriate wells of the Glutamate Microtiter Strips .
2.	Pipette 50 µl of the Glutamate Antiserum into all wells and mix shortly.
3.	Cover plate with Adhesive Foil and incubate for 15 – 20 h (overnight) at 2 – 8 °C .
4.	Remove the foil. Discard or aspirate the content of the wells. Wash the plate 3 x by adding 300 µl of Wash Buffer , discarding the content and blotting dry each time by tapping the inverted plate on absorbent material.
5.	Pipette 100 µl of the Enzyme Conjugate into all wells.
6.	Incubate for 30 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm).
7.	Discard or aspirate the contents of the wells and wash the plate 3 x by adding 300 µl of Wash Buffer , discarding the content and blotting dry each time by tapping the inverted plate on absorbent material.
8.	Pipette 100 µl of the Substrate into all wells and incubate for 20 – 30 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). Avoid exposure to direct sunlight!
9.	Add 100 µl of the Stop Solution to each well and shake the microtiter plate to ensure a homogeneous distribution of the solution.
10.	Read the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to 450 nm (if available a reference wavelength between 620 nm and 650 nm is recommended).

7. Calculation of results

Measuring range	Glutamate	
	Urine	0.26 – 60 µg/ml

The standard curve is obtained by plotting the absorbance readings (calculate the mean absorbance) of the standards (linear, y-axis) against the corresponding standard concentrations (logarithmic, x-axis). Use non-linear regression for curve fitting (e. g. 4-parameter, marquardt).

⚠ *This assay is a competitive assay. This means: the OD-values are decreasing with increasing concentrations of the analyte. OD-values found below the standard curve correspond to high concentrations of the analyte in the sample and have to be reported as being positive.*

Urine samples and controls

The concentrations of the samples and controls can be read directly from the standard curve. Samples found with concentrations higher than the highest standard (Standard F) should be diluted accordingly with water (deionized, distilled, or ultra-pure) and have to be re-assayed.

Conversion

Glutamate ($\mu\text{g/ml}$) $\times 6.8 =$ Glutamate ($\mu\text{mol/l}$)

Expected reference values

It is strongly recommended that each laboratory should determine its own reference values.

Spontaneous urine
1034 – 7726 $\mu\text{g/g}$ creatinine
7 – 52.5 $\mu\text{mol/g}$ creatinine
0.8 – 5.9 mmol/mol creatinine


Values significantly outside the reference range should be assessed by a doctor.

7.1 Quality control

The confidence limits of the kit controls are indicated on the QC-Report.

7.2 Typical standard curve

⚠ Example, do not use for calculation!

8. Assay characteristics

Analytical Sensitivity	Glutamate
Limit of Blank (LOB)	0.11 $\mu\text{g/ml}$
Limit of Detection (LOD)	0.17 $\mu\text{g/ml}$
Limit of Quantification (LOQ)	0.26 $\mu\text{g/ml}$

Analytical Specificity (Cross Reactivity)	Substance	Cross Reactivity (%)
		Glutamate
	L-Glutamine	< 0.4
	Glycine	< 0.4
	β -Alanine	< 0.4
	L-Alanine	< 0.4
	L-Aspartic Acid	< 0.4
	GABA	< 0.4
	5-Amino-n-valeric Acid	< 0.4

Precision							
Intra-Assay				Inter-Assay			
Sample	n	Mean $\pm$ SD ($\mu\text{g/ml}$)	CV (%)	Sample	n	Mean $\pm$ SD ($\mu\text{g/ml}$)	CV (%)
1	10	0.8 ± 0.1	10.8	1	13	1.7 ± 0.24	14.3
2	10	1.3 ± 0.1	8.7	2	14	5.0 ± 0.57	11.4
3	10	2.2 ± 0.1	6.3	3	14	10.6 ± 0.73	6.9
4	10	4.8 ± 0.2	4.0	4	13	3.0 ± 0.43	14.2
5	10	12.5 ± 0.6	4.6	5	14	5.6 ± 0.71	12.5
6	10	39.7 ± 2.2	5.6	6	14	10.0 ± 0.87	8.7

Linearity		Serial dilution up to	Range (%)	Mean (%)
	Urine	1:64	94 – 113	105

Recovery		Range ($\mu\text{g/ml}$)	Range (%)	Mean (%)
	Urine	1.25 – 41.0	97 – 108	102

9. References/Literature

- (1) Bieger, W.P., NeuroStress Guide. 2011.
- (2) Bustillo, J.R., Use of proton magnetic resonance spectroscopy in the treatment of psychiatric disorders: a critical update. Dialogues Clin Neurosci, 2013. 15(3): p. 329-37.
- (3) Duman, R.S., G. Sanacora, and J.H. Krystal, Altered Connectivity in Depression: GABA and Glutamate Neurotransmitter Deficits and Reversal by Novel Treatments. Neuron, 2019. 102(1): p. 75-90.
- (4) Femenia, T., et al., Dysfunctional hippocampal activity affects emotion and cognition in mood disorders. Brain Res, 2012. 1476: p. 58-70.
- (5) Flasnoecker, M., Reise aus dem Stress - Körper, Geist und Psyche stärken. 2015: W. Zuckschwerdt Verlag GmbH.
- (6) Frisardi, V., F. Panza, and A.A. Farooqui, Late-life depression and Alzheimer's disease: the glutamatergic system inside of this mirror relationship. Brain Res Rev, 2011. 67(1-2): p. 344-55.
- (7) Gao, S.F. and A.M. Bao, Corticotropin-releasing hormone, glutamate, and gamma-aminobutyric acid in depression. Neuroscientist, 2011. 17(1): p. 124-44.
- (8) Harris, R.E. and D.J. Clauw, Imaging central neurochemical alterations in chronic pain with proton magnetic resonance spectroscopy. Neurosci Lett, 2012. 520(2): p. 192-6.
- (9) Kendell, S.F., J.H. Krystal, and G. Sanacora, GABA and glutamate systems as therapeutic targets in depression and mood disorders. Expert Opin Ther Targets, 2005. 9(1): p. 153-68.
- (10) Krystal, J.H., et al., Glutamate and GABA systems as targets for novel antidepressant and mood-stabilizing treatments. Mol Psychiatry, 2002. 7 Suppl 1: p. S71-80.
- (11) Lener, M.S., et al., Glutamate and Gamma-Aminobutyric Acid Systems in the Pathophysiology of Major Depression and Antidepressant Response to Ketamine. Biol Psychiatry, 2017. 81(10): p. 886-897.
- (12) Sanacora, G., G. Treccani, and M. Popoli, Towards a glutamate hypothesis of depression: an emerging frontier of neuropsychopharmacology for mood disorders. Neuropharmacology, 2012. 62(1): p. 63-77.
- (13) Strienz, J., Nebennierenunterfunktion - Stress stört die Hormonbalance. 3 ed. 2019: W. Zuckschwerdtverlag München.

For updated literature or any other information please contact your local supplier.

Symbols:

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Use-by date		Batch code		For in-vitro diagnostic use only!
	Consult instructions for use		Content		CE marking of conformity
	Caution		Catalogue number		Distributor
	Date of manufacture				

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9. Referenzen/Literatur

- (1) Bieger, W.P., NeuroStress Guide. 2011.
- (2) Bustillo, J.R., Use of proton magnetic resonance spectroscopy in the treatment of psychiatric disorders: a critical update. Dialogues Clin Neurosci, 2013. 15(3): p. 329-37.
- (3) Duman, R.S., G. Sanacora, and J.H. Krystal, Altered Connectivity in Depression: GABA and Glutamate Neurotransmitter Deficits and Reversal by Novel Treatments. Neuron, 2019. 102(1): p. 75-90.
- (4) Femenia, T., et al., Dysfunctional hippocampal activity affects emotion and cognition in mood disorders. Brain Res, 2012. 1476: p. 58-70.
- (5) Flasnoecker, M., Reise aus dem Stress - Körper, Geist und Psyche stärken. 2015: W. Zuckschwerdt Verlag GmbH.
- (6) Frisardi, V., F. Panza, and A.A. Farooqui, Late-life depression and Alzheimer's disease: the glutamatergic system inside of this mirror relationship. Brain Res Rev, 2011. 67(1-2): p. 344-55.
- (7) Gao, S.F. and A.M. Bao, Corticotropin-releasing hormone, glutamate, and gamma-aminobutyric acid in depression. Neuroscientist, 2011. 17(1): p. 124-44.
- (8) Harris, R.E. and D.J. Clauw, Imaging central neurochemical alterations in chronic pain with proton magnetic resonance spectroscopy. Neurosci Lett, 2012. 520(2): p. 192-6.
- (9) Kendell, S.F., J.H. Krystal, and G. Sanacora, GABA and glutamate systems as therapeutic targets in depression and mood disorders. Expert Opin Ther Targets, 2005. 9(1): p. 153-68.
- (10) Krystal, J.H., et al., Glutamate and GABA systems as targets for novel antidepressant and mood-stabilizing treatments. Mol Psychiatry, 2002. 7 Suppl 1: p. S71-80.
- (11) Lener, M.S., et al., Glutamate and Gamma-Aminobutyric Acid Systems in the Pathophysiology of Major Depression and Antidepressant Response to Ketamine. Biol Psychiatry, 2017. 81(10): p. 886-897.
- (12) Sanacora, G., G. Treccani, and M. Popoli, Towards a glutamate hypothesis of depression: an emerging frontier of neuropsychopharmacology for mood disorders. Neuropharmacology, 2012. 62(1): p. 63-77.
- (13) Strienz, J., Nebennierenunterfunktion - Stress stört die Hormonbalance. 3. ed. 2019: W. Zuckschwerdtverlag München.

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